# **Table of Contents**

#### **General Information**

Gigahertz-Optik manufactures light detectors for use in combination with Gigahertz-Optik optometers and accessories to form complete radiometric and photometric measurement systems. Detectors can also be used with integrating spheres to form the basis of uniform light sources and reflectance/ transmittance instruments. A typical sensor head may consist of a photodiode, optical bandpass filter, diffuser, lens assembly or other optical elements. A complete range of detector head housings and configurations designed around the specific application demands of dynamic range, spectral range, physical size and measurement environment are available. All optical radiation detectors are calibrated and certified by Gigahertz-Optik's calibration laboratory with both absolute sensitivity data and spectral responsivity plot provided.

Table of Cont	ents								
Light Detecto	r Connector Types for D	ifferent G.O. Op	otometers						Page 62
Calculation of	f the Maximum Measura	ble Intensity an	d Maximum	Attainable I	Resolution				Page 63
Model	Spectral Function	Wavelength Range	Diffuser Diameter	Aperture (mm)	Sphere Diameter	Unit	Pac	kage	Specifi- cations
		nm	mm	F.O.V.	mm		Туре	Page	Page
		Com	npact Size II	luminance [	Detectors				
VL-1101	Photopic; f₁≤5%	380-780	-	1,1 mm Ø	-	lx, fc	11	90	63
		Photopi	ic & Scotop	ic Illuminan	ce Detectors				
VL-3701	Photopic; f₁≤3%	380-780	7	Cosine	-	lx, fc	37/20	90	64
VL-3702	Photopic; f₁≤5%	380-780	7	Cosine	-	lx, fc	37/20	90	64
VL-3704	Photopic; f₁≤5%	380-780	7	Cosine	-	lx, fc	37/20	90	64
VL-3705	Scotopic; f₁≤5%	380-780	7	Cosine	-	lx, fc	37/20	90	64
PD-9310	Photopic; f₁≤3 or 6%	380-780	30	Cosine	-	lx, fc	PD93	91	64
		Temperat	ture Control	lled Illumina	nce Detector	rs			
TD-11VL01	Photopic; f₁≤3%	380-780	-	-	-	lx, fc	TD-11	90	65
	- <b>·</b>	Compact	Size Broad	band Irradia	nce Detecto	rs			
RW-1103	Radiometric VISIBLE	400-800	-	1.1 x 1.1	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	11	90	66
RW-1104	Radiometric NIR	800-1000	-	1.1 x 1.1	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	11	90	66
RW-1105	Radiometric VISNIR	400-1000	-	1.1 x 1.1	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	11	90	66
	- <b>·</b>	Temperature C	Controlled B	roadband Ir	radiance De	tectors			
TD-11RW03	Radiometric VISIBLE	400-800	-	2.8 x 2.8	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	TD-11	90	66
TD-11RW04	Radiometric NIR	800-1000	-	2.8 x 2.8	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	TD-11	90	66
TD-11RW05	Radiometric VISNIR	400-1000	-	2.8 x 2.8	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	TD-11	90	66
		High Sens	sitive Broad	band Irradia	ance Detecto	rs			
RW-3701	Radiometric BLUE	400-500	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	67
RW-3702	Radiometric RED	700-800	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	67
RW-3703	Radiometric VISIBLE	400-800	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	67
RW-3704	Radiometric NIR	800-1000	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	67
RW-3705	Radiometric VISNIR	400-1000	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	67
	- <b>-</b>	Compact S	ize Broadba	and UV Irrad	iance Detect	ors			
UV-1101	Radiometric UV-A	315-400	-	2.5 x 2.5	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	11	90	68
UV-1102	Radiometric UV-B	280-315	-	1.1 x 1.1	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	11	90	68
	Т	emperature Co	ntrolled Bro	adband UV	Irradiance D	etectors			
TD-11UV01	Radiometric UV-A	315-400	-	2.5 x 2.5	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	TD-11	90	68
TD-11UV02	Radiometric UV-B	280-315	-	1.1 x 1.1	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	TD-11	90	68
		High Sensiti	ivity Broadb	and UV Irra	diance Detec	tors			
UV-3701	Radiometric UV-A	315-400	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69
UV-3702	Radiometric UV-B	280-315	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/50	90	69
UV-3703	Radiometric UV-C	250-280	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69

# Page 60 - Light Detectors



# Table of Contents

Model	Spectral Function	Wavelength Range	Diffuser Diameter	Aperture (mm)	Sphere diameter	Unit	Pack	age	Specifi- cations
		nm	mm	F.O.V.	mm		Туре	Page	Page
		High Sensitivit	y Broadban	d UV Irradia	ance Detect	ors			
UV-3710	Radiometric UV-A	320-400	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69
UV-3711	Radiometric UV-B	280-320	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/50	90	69
UV-3716	Radiometric UV-A	305-400	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69
UV-3717	UV-A Visible Blind	325-400	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69
UV-3718	Narrowband UV-C	254	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69
UV-3719	Broadband UV	250-400	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69
		Photobiolo	gical Actini	c Irradiance	Detectors		1		1
UV-3704	UV-Erythema	250-400	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69/70
UV-3706	Bilirubin	420-530	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69/70
UV-3708	ACGIH/ICNIRP	250-400	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/50	90	69/70
UV-3709	Blue-Light Hazard	400-520	11	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/32	90	69/70
	<u> </u>	Plant Pl	hysiology Ir	radiance De	etectors				
PS-3701	PAR	400-700	11	Cosine	-	µMol/m²/s	37/32	90	71
PS-3702	Photon Effective	320-500	11	Cosine	-	μMol/m²/s	37/32	90	71
PS-3703	Photon Effective	590-900	11	Cosine	-	μMol/m²/s	37/32	90	71
	High Int	ensity/Temperat	ure Irradian	ce Detecto	rs / UV-Curi	1	,		1
RCH-006	Broadband UV	250-400	9	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/8/F	91	72/73
RCH-008	Narrowband UVA	365 peak	9	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/8/F	91	72/73
RCH-009	Narrowband BLUE	430 peak	9	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/8/F	91	72/73
RCH-106	Broadband UV	250-400	9	Cosine	_	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/8/R	91	72/73
RCH-108	Narrowband UVA	365 peak	9	Cosine	_	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/8/R	91	72/73
RCH-109	Narrowband BLUE	430 peak	9	Cosine	_	W/m <sup>2</sup> ; W/cm <sup>2</sup>	37/8/R	91	72/73
RCH-506	Broadband UV	250-400	6	Cosine	_	W/m <sup>2</sup> ; W/cm <sup>2</sup>	15	91	72/73
RCH-508	Narrowband UVA	365 peak	6	Cosine	_	W/m <sup>2</sup> ; W/cm <sup>2</sup>	15	91	72/73
RCH-509	Narrowband BLUE	430 peak	6	Cosine		W/m <sup>2</sup> ; W/cm <sup>2</sup>	15	91	72/73
	Nanowband BECE	360 Degree Fie			ance Detect		10	51	12,10
ROD-360	Narrowband UV-C	254	ROD	360 Deg		W/m <sup>2</sup> : W/cm <sup>2</sup>	ROD	92	73
NOD-360	Narrowband 0V-C	Multi-Channe		0	- Dotocto	, , , ,	nob	92	73
XD-9501	UV-A 315-400 & UV		8	Cosine		[	XD95	91	74/75
XD-9501 XD-9502	UVA 315-400 & 0V		8	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup> W/m <sup>2</sup> ; W/cm <sup>2</sup> & lx	XD95	91	74/75 74/75
XD-9502 XD-9503		•			-		XD95		
XD-9503 XD-9504	UV-A/B 240-340 & U Radiometric VISIE		8	Cosine Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup> W/m <sup>2</sup> ; W/cm <sup>2</sup>	XD95	91 91	74/75 74/75
XD-9505	UV Broadband	250-400		Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	XD95	91	74/75
XD-9505 XD-9506			8		-				
XD-9506 XD-9509	ACGIH 250-314/315-400		8	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup> W/m <sup>2</sup> ; W/cm <sup>2</sup>	XD95	91	74/75
	Erythema I250-32		8	Cosine	-		XD95	91	74/75
XD-9510	ACGIH 250-314/315-400		8	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	XD95	91	74/75
XD-9511	UV-A visible blind	315-400	8	Cosine	-	W/m <sup>2</sup> ; W/cm <sup>2</sup>	XD95	91	74/75
TP-4501	Triple PAR & Pl	•	8	Cosine	-	µMol/m²/s & lx	45	91	74/75
			s Color & III	1	Jetectors		0		
CT-3701	x,y,z, 2 Degrees	380-780	10	Cosine	-	x,y; u,v; K, lx	37/TC	92	76
CT-4501	x,y,z, 2 Degrees	380-780	10	Cosine	-	x,y; u,v; K, lx	45	90	76
			nance & Rad						
LDM-9810	Interchangeable PD- Spectral function PD-			° and 6° Sel		PD-16 dependent	LDM98	92	77/78
LDM-9811			1.7, 11 ar	nd 100 mrad	Selectable	PD-16 dependent	LDM98	92	77/78
PD-16-BLH	Blue Light Hazard	400-520	-	5.8 x 5.8	-	W/m <sup>2</sup> sr	PD16	92	77/78

# Light Detectors - Page 61

Table of Contents

LSM-9901     Photopic,       PRW-0505     Radiometric       Integrating     Spheres from       Spheres     Si Photocic       LP-9901     Si Photocic       LP-0101     InGaAs Pho       LP-0102     Si & InGaAs Pho       LP-0103     Si Photocic       LP-9901     Si Photocic       LP-9910     Si Photocic       LP-9920     InGaAs Pho       LP-9930     Si & InGaAs Pho       LP-9930     Si & InGaAs Pho       LP-9930     Si & InGaAs Pho       PD-1101     Si Photocic       PD-1102     InGaAs Pho       PD-1103     InGaAs Pho       PD-1104     InGaAs Pho       PD-1105     Si LP Photo       PD-1106     Si LP Photo       PD-1107     Si LP Photo       PD-1108     Si LP Photo       PD-1109     GaP Photo       PD-1108     Si LP Photo       PD-1109     GaP Photo       PD-1101     Si Photoc       PD-1102     InGaAs Pho       PD-1103     Si LP Photoc	I Function	Wavelength Range	Diffuser Diameter	Aperture (mm)	Sphere Diameter	Unit	Pack	age	Specifi- cation
PD-16-RTHARTH IFPD-16-RW05RadiometricPD-16-VL01Photopic,LDM-98Z-NL-LDM-98Z-FC-LDM-9901Photopic,SRT-M37A Range ofUSM-9901Photopic,RW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si Photopic,LP-9901Si Photopic,PRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si Photopic,LP-9901Si Photopic,LP-0102Si & InGaAs Pho Si LP Photopic,LP-0103Si Photopic,LP-0104InGaAs Pho Si LP Photopic,LP-9930Si & InGaAs Pho Si LP Photopic,PD-1101Si Photopic,PD-1102InGaAs Pho Photopic,PD-1103InGaAs Pho Si LP Photopic,PD-1104InGaAs Pho Si LP Photopic,PD-1105Si LP Photopic,PD-1106Si LP Photopic,PD-1107Si LP Photopic,PD-1108Si LP Photopic,PD-1109GaP Photopic,PD-1109GaP Photopic,PD-1101Si Photopic,PD-1103InGaAs Pho Si LP Photopic,PD-1104InGaAs Photopic,PD-1105Si LP Photopic,PD-1106Si LP Photopic,PD-1107Si LP Photopic,PD-1108Si LP Photopic,PD-1109GaAs Photopic,PD-1101Si Photopic,PD-1103InGaAs Photopic,PD		nm	mm	F.O.V.	mm		Туре	Page	Page
PD-16-RTHARTH IFPD-16-RW05RadiometricPD-16-VL01Photopic,LDM-98Z-NL-LDM-98Z-FC-LDM-9901Photopic,SRT-M37A Range ofUSM-9901Photopic,RW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si Photopic,LP-9901Si Photopic,PRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si Photopic,LP-9901Si Photopic,LP-0102Si & InGaAs Pho Si PhotocLP-0103Si PhotocLP-9910Si LP PhotoLP-9920InGaAs Pho Si LP PhotoLP-9930Si & InGaAs Pho Photopic,PD-1101Si PhotocPD-1102InGaAs Pho Photopic,PD-1103InGaAs Pho Si LP PhotocPD-1104InGaAs Pho Si LP PhotocPD-1105Si LP Photoc Si LP PhotocPD-1106Si LP Photoc Si LP PhotocPD-1107Si LP Photoc Si LP PhotocPD-1108Si LP Photoc Si LP PhotocPD-1109GaAs Pho Si LP PhotocPD-1109Spheres from SpheresTD-1101Si Photoc Si PhotocTD-1102InGaAs Pho Spheres from SpheresTD-1103InGaAs Pho Spheres from SpheresTD-1103Si Photoc Spheres from Spheres from Spheres from		Lur	ninance & R	adiance Dete	ectors		1		
PD-16-RW05RadiometricPD-16-VL01Photopic, ?LDM-98Z-FC-LDM-9901Photopic, ?SRT-M37A Range ofUN-9901Photopic, ?SRT-M37A Range ofUN-9901Photopic, ?RW-0505RadiometricIntegrating SpheresSpheres from RWUP-0101InGaAs PhotopicLP-0101Si & InGaAs PhotopicLP-0102Si & InGaAs PhotopicLP-0103Si PhotocicLP-0103Si PhotocicLP-9920InGaAs PhotopicLP-9930Si & InGaAs PhotopicLP-9910Si LP PhotopicPD-1101Si LP PhotopicPD-1102InGaAs PhotopicPD-1103Si LP PhotopicPD-1104InGaAs PhotopicPD-1105Si LP PhotopicPD-1106Si LP PhotopicPD-1107Si LP PhotopicPD-1108Si LP PhotopicPD-1109GaP PhotopicPD-1109Si LP PhotopicPD-1101Si LP PhotopicPD-1103Si LP PhotopicPD-1104Si LP PhotopicPD-1105Si LP PhotopicPD-1106Si LP PhotopicPD-1107Si LP PhotopicPD-1108Si LP PhotopicPD-1109GaAsPhopicPD-1103InGaAsPhopicPD-1104Si LP PhotopicPD-1105Si LP PhotopicPD-1105Si LP PhotopicPD-1105Si LP PhotopicPD-1105Si LP PhotopicPD-1105 <t< td=""><td>ermal Hazard</td><td>500-1200</td><td>-</td><td>5.8 x 5.8</td><td>-</td><td>W/m²sr</td><td>PD16</td><td>92</td><td>77/78</td></t<>	ermal Hazard	500-1200	-	5.8 x 5.8	-	W/m²sr	PD16	92	77/78
PD-16-VL01Photopic, iLDM-98Z-NL-LDM-98Z-FCPhotopic, iSRT-M37A Range ofLSM-9901Photopic, iPRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si PhotooicLP-0102Si & InGaAs Photopic, iLP-0103Si PhotooicLP-0103Si PhotooicLP-9901Si A InGaAs Photopic, iLP-0103Si PhotooicLP-9910Si PhotooicLP-9920InGaAs Photopic, iLP-9930Si & InGaAs Photopic, iLP-9930Si & InGaAs Photopic, iLP-9930Si & InGaAs Photopic, iLP-9910Si LP Photopic, iPD-1103InGaAs Photopic, iPD-1104InGaAs Photopic, iPD-1105Si LP Photopic, iPD-1106Si LP Photopic, iPD-1107Si LP Photopic, iPD-1108Si LP Photopic, iPD-1109GaP Photopic, iPD-1109GaAs Photopic, iPD-1101Si Photopic, iPD-1102InGaAs Photopic, iPD-1103Si LP Photopic, iPD-1104Si LP Photopic, iPD-1105Si LP Photopic, iPD-1106Si LP Photopic, iPD-1107Si LP Photopic, iPD-1108Si LP Photopic, iPD-1109GaAs Photopic, iPD-1101Si Photopic, iPD-1102InGaAs Photopic, iPD-1103Si Photopic, iPD-1104Si LP Photopic, iPD-1105	H IR-A	800-1200	-	5.8 x 5.8	-	W/m²sr	PD16	92	77/78
LDM-98Z-NL-LDM-98Z-FC-LDM-9901Photopic,SRT-M37A Range ofLSM-9901Photopic,PRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si Photopic,LP-9901Si Photopic,LP-0102Si A InGaAs Pho LP-0103LP-0103Si PhotocLP-9901Si LP Photo Si LP PhotoLP-9901Si LP Photo Si LP PhotoLP-9910Si A InGaAs Pho LP-9920LP-9920InGaAs Pho Si LP PhotoLP-9930Si & InGaAs Pho D-1102PD-1101Si Photoc Si LP PhotoPD-1102InGaAs Pho D Photopic,PD-1103InGaAs Pho Photopic,PD-1104Si LP Photopic,PD-1105Si C Photopic,PD-1106Si LP Photopic,PD-1107Si LP Photopic,PD-1108Si LP Photopic,PD-1109GaAs Pho CPhotopic,PD-1109Si A DASP Pho Si LP Photopic,PD-1101Si LP Photopic,PD-1103InGaAs Pho Si LP Photopic,PD-1104Si LP Photopic,PD-1105Si LP Photopic,PD-1106Si LP Photopic,PD-1107Si LP Photopic,PD-1108Si LP Photopic,PD-1109GaAs Pho Si LP Photopic,PD-1103InGaAs Pho Si LP Photopic,PD-1104Si LP Photopic,PD-1105Si LP Photopic,PD-1105Si LP Photopic,PD-1105Si LP Photop	etric VIS-NIR	400-1000	-	5.8 x 5.8	-	W/m²sr	PD16	92	77/78
LDM-98Z-FC-LDM-9901Photopic, "SRT-M37A Range ofUSM-9901Photopic, "PRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si PhotocLP-0102Si & InGaAs PhotocLP-0103Si PhotocLP-9901Si PhotocLP-9901Si PhotocLP-9910Si PhotocLP-9920InGaAs PhotocLP-9930Si & InGaAs PhotocLP-9910Si LP PhotocPD-1101Si PhotocPD-1102InGaAs PhoPD-1103Si LP PhotocPD-1104Si LP PhotocPD-1105Si LP PhotocPD-1106Si LP PhotocPD-1107Si LP PhotocPD-1108Si LP PhotocPD-1109GaP PhotocPD-1109GaP PhotocPD-1101Si PhotocTD-1102InGaAs PhoIntegrating Spheres fromIntegratingSpheres fromIntegratingSpheres fromIntegratingSpheres from	oic, f1≤3%	380-780	-	5.8 x 5.8	-	cd/m <sup>2</sup>	PD16	92	77/78
LDM-9901Photopic,SRT-M37A Range ofLSM-9901Photopic,PRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si PhotoolLP-0102Si & InGaAs Photopic,LP-0103Si PhotoolLP-0103Si & InGaAs Photopic,LP-9910Si PhotoolLP-9920InGaAs Photopic,LP-9920InGaAs Photopic,LP-9930Si & InGaAs Photopic,LP-9910Si PhotoolLP-9920InGaAs Photopic,LP-9930Si & InGaAs Photopic,LP-1101Si LP Photopic,PD-1102InGaAs Photopic,PD-1103InGaAs Photopic,PD-1104InGaAs Photopic,PD-1105Si LP Photopic,PD-1106Si LP Photopic,PD-1107Si LP Photopic,PD-1108Si LP Photopic,PD-1109GaP Photopic,PD-1109GaP Photopic,PD-1101Si Photopic,PD-1103Si LP Photopic,PD-1104InGaAs Photopic,PD-1105Si LP Photopic,PD-1106Si LP Photopic,PD-1107Si LP Photopic,PD-1108Si LP Photopic,PD-1109GaAsP Photopic,PD-1109Spheres fromSpheresSpheres fromSpheresSi Photopic,TD-1102InGaAs Photopic,TD-1103InGaAs Photopic,TD-1104Si Photopic,SpheresSpheres fromSpheresSpheres from	-	380 - 1000	-	50	-	-	LDM98	92	77/78
SRT-M37A Range ofLSM-9901Photopic, iPRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si PhotocicLP-9901Si PhotocicLP-0102Si & InGaAs PhoLP-0103Si PhotocicLP-0103Si PhotocicLP-9901Si & InGaAs PhoLP-9910Si PhotocicLP-9920InGaAs PhoLP-9930Si & InGaAs PhoLP-9930Si & InGaAs PhoLP-9930Si & InGaAs PhoPD-1101Si PhotocicPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotocicPD-1106Si LP PhotocicPD-1107Si LP PhotocicPD-1108Si LP PhotocicPD-1109GaP PhotocicPD-1109GaP PhotocicPD-1109Spheres fromSpheresSpheres fromSpheresSpheres fromSpheresSpheres fromSpheresSpheres fromSpheresSpheres fromSpheresSpheres fromSpheresSpheres fromSpheresSpheres from	-	Fiber Coupler f	or flexible ligh	nt guide bundl	es			92	77/78
LSM-9901 Photopic, PRW-0505 Radiometric Integrating Spheres from Spheres Si Photoco LP-0101 Si Photoco LP-0102 Si & InGaAs Pho LP-0103 Si Photoco LP-0201 Si LP Photoco LP-9910 Si Photoco LP-9920 InGaAs Pho LP-9920 Si & InGaAs Pho LP-9930 Si & InGaAs Pho LP-9930 Si & InGaAs Pho PD-1101 Si LP Photoco PD-1102 InGaAs Pho PD-1102 InGaAs Pho PD-1103 InGaAs Pho PD-1104 InGaAs Pho PD-1105 SiC Photoco PD-1105 Si LP Photoco PD-1106 Si LP Photoco PD-1107 Si LP Photoco PD-1108 Si LP Photoco PD-1109 GaP Photoco PD-1109 GaP Photoco PD-1109 GaP Photoco PD-1109 Spheres from Spheres from Spheres from TD-1103 InGaAs Pho	oic, f1≤5%	380-780	-	1°	-	cd/m <sup>2</sup>	LDM99	92	79
PRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si PhotocLP-9901Si PhotocLP-0102Si & InGaAs Pho Si & InGaAs Pho LP-0103Si PhotocLP-0103Si PhotocLP-0103Si PhotocLP-9910Si PhotocLP-9920InGaAs Pho Si & InGaAs Pho LP-9930LP-9930Si & InGaAs Pho PhotocPD-1101Si PhotocPD-1102InGaAs Pho PhotocPD-1103InGaAs Pho PhotocPD-1104InGaAs Pho PhotocPD-1105SiC Photoc PhotocPD-1106Si LP Photoc PhotocPD-1107Si LP Photoc Si C PhotocPD-1108Si LP Photoc PhotocPD-1109GaAs Pho SpheresTD-1101Si Photoc Spheres from SpheresTD-1102InGaAs Pho Spheres from Spheres fromTD-1103InGaAs Pho Spheres from Spheres from	e of Front Lens	ses is Offered fo	r the VL, RW	and UV Detect	tors for Lumir	nance & Radiance S	et-ups	89	80
PRW-0505RadiometricIntegrating SpheresSpheres from RWLP-9901Si PhotocLP-9901Si PhotocLP-0102Si & InGaAs Pho LP-0103LP-0103Si PhotocLP-0103Si PhotocLP-9910Si LP PhotocLP-9920InGaAs Pho LP-9930LP-9930Si & InGaAs Pho LP-9930LP-9930Si & InGaAs Pho PD-1101PD-1101Si PhotocPD-1102InGaAs Pho PD-1103PD-1103InGaAs Pho PD-1104PD-1104InGaAs Pho PD-1105PD-1105Si LP Photoc PD-1107PD-1108Si LP Photoc PD-1108PD-1109GaP Photoc PD-1109PD-1101Si Photoc PD-1105PD-1103InGaAs Pho Photoc PD-1109PD-1104Si LP Photoc Photoc PD-1105PD-1105Si LP Photoc Photoc PD-1109PD-1108Si LP Photoc Photoc PD-1109PD-1109GaAsP Photoc Photoc PD-1101Integrating TD-1102Spheres from 		Lumino	us Flux & Ra	adiant Power	Detectors				
Integrating SpheresSpheres from RWLP-9901Si PhotodLP-0101InGaAs PhoLP-0102Si & InGaAs PhoLP-0103Si PhotodLP-0103Si PhotodLP-9910Si PhotodLP-9920InGaAs PhoLP-9930Si & InGaAs PhoLP-9910Si LP PhotodLP-9920InGaAs PhoLP-9930Si & InGaAs PhoPD-1101Si PhotodPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotodPD-1106Si LP PhotodPD-1107Si LP PhotodPD-1108Si LP PhotodPD-1109GaP PhotodPD-1109GaP PhotodPD-1109Spheres fromSpheresSpheres fromTD-1101Si PhotodTD-1103InGaAs PhoTD-1103Spheres from	 bic, f1≤5%	380-780	-	12.5 Ø	50	lm	11/50	91	81
Spheres     RW       LP-9901     Si Photos       LP-0101     InGaAs Pho       LP-0102     Si & InGaAs Pho       LP-0103     Si Photos       LP-0103     Si Photos       LP-0103     Si Photos       LP-0201     Si LP Photos       LP-9910     Si Photos       LP-9920     InGaAs Pho       LP-9930     Si & InGaAs Pho       LP-9930     Si & InGaAs Pho       PD-1101     Si Photos       PD-1102     InGaAs Pho       PD-1103     InGaAs Pho       PD-1104     InGaAs Pho       PD-1105     Si LP Photos       PD-1106     Si LP Photos       PD-1107     Si LP Photos       PD-1108     Si LP Photos       PD-1109     GaP Photos       PD-1109     GaP Photos       PD-1109     GaP Photos       PD-1109     GaP Photos       PD-1101     Si Photos       PD-1102     InGaAs Pho       Integrating     Spheres from       Spheres     InGaAs Pho	etric VIS-NIR	400-1000	-	12.5 Ø	50	W	11/50	91	81
LP-9901     Si Photoco       LP-0101     InGaAs Pho       LP-0102     Si & InGaAs Pho       LP-0103     Si Photoco       LP-0201     Si LP Photoco       LP-9910     Si Photoco       LP-9920     InGaAs Pho       LP-9930     Si & InGaAs Pho       LP-9930     Si & InGaAs Pho       LP-9930     Si & InGaAs Pho       PD-1101     Si LP Photoco       PD-1102     InGaAs Pho       PD-1103     InGaAs Pho       PD-1104     InGaAs Pho       PD-1105     Si LP Photoco       PD-1106     Si LP Photoco       PD-1107     Si LP Photoco       PD-1108     Si LP Photoco       PD-1109     GaP Photoco       PD-1108     Si LP Photoco       PD-1109     GaP Photoco       PD-1109     GaP Photoco       PD-1109     GaP Photoco       PD-1101     Si Photoco       PD-1102     InGaAs Pho       PD-1103     Si Photoco       PD-1104     GaAs Pho       Integrating     Spheres	from 100 mm	(4 inch) to 500 n	nm (20 inch)	Diameter are (	Dffered which	can be Combined	with VL,	See Ir	ntegrating
LP-0101InGaAs PhoLP-0102Si & InGaAs PhoLP-0103Si PhotocLP-0103Si PhotocLP-0201Si LP PhotoLP-9910Si PhotocLP-9920InGaAs PhoLP-9930Si & InGaAs PhoLP-9930Si & InGaAs PhoLP-9910Si LP PhotoLP-9920InGaAs PhoPD-1101Si LP PhotoPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhoIntegrating SpheresSpheres fromTD-1103InGaAs PhoInGaAs PhoInGaAs PhoSpheresSpheres from						surement Set-ups		Sp	oheres
LP-0101InGaAs PhoLP-0102Si & InGaAs PhoLP-0103Si PhotocLP-0103Si PhotocLP-0201Si LP PhotoLP-9910Si PhotocLP-9920InGaAs PhoLP-9930Si & InGaAs PhoLP-9930Si & InGaAs PhoLP-9910Si LP PhotoLP-9920InGaAs PhoPD-1101Si LP PhotoPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhoIntegrating SpheresSpheres fromTD-1103InGaAs PhoInGaAs PhoInGaAs PhoSpheresSpheres from		Laser St	ray-light and	d Laser Powe	r Detectors				
LP-0102Si & InGaAs PLP-0103Si PhotocLP-0201Si LP PhotoLP-9910Si PhotocLP-9920InGaAs PhoLP-9930Si & InGaAs PLP-9930Si & InGaAs PhoLP-9930Si & InGaAs PhoPD-1101Si PhotocPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105Si LP PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109Spheres fromSpheresSpheres fromTD-1102InGaAs PhoTD-1103InGaAs Pho	otodiode	350-1100	-	7 Ø	-	W/m² nm / W nm	37/8.5	91	82
LP-0102Si & InGaAs PLP-0103Si PhotocLP-0201Si LP PhotoLP-9910Si PhotocLP-9920InGaAs PhoLP-9930Si & InGaAs PLP-9930Si & InGaAs PhoLP-9930Si & InGaAs PhoPD-1101Si PhotocPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105Si C PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109Spheres fromSpheresSpheres fromTD-1102InGaAs PhoTD-1103InGaAs PhoTD-1103Spheres from		Laser Pow	er Detectors	s with Integra	ting Sphere	s			
LP-0103Si PhotoLP-0201Si LP PhotoLP-9910Si PhotoLP-9920InGaAs PhoLP-9930Si & InGaAs PhoLP-9930Si & InGaAs PhoLP-9930Si & InGaAs PhoPD-1101Si LP PhotoPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109Spheres fromSpheresSpheres fromTD-1103InGaAs PhoTD-1103InGaAs PhoTD-1103Spheres from	Photodiode	800-1800	-	5 Ø	30	W nm	11/30	91	82/83
LP-0201Si LP PhotoLP-9910Si PhotocoLP-9920InGaAs PhoLP-9930Si & InGaAs PhoLPPA-9901Si LP PhotoPD-1101Si PhotocoPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105Si LP PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109Spheres fromSpheresSpheres fromTD-1103InGaAs PhoInGaAs PhoSi LP PhotoSpheresSpheres fromSpheresSpheres fromTD-1103InGaAs PhoIntegratingSpheres from	s Photodiode	350-1800	-	5 Ø	30	W nm	11/30	91	82/83
LP-9910Si PhotocLP-9920InGaAs PhoLP-9930Si & InGaAs PhoLP9930Si & InGaAs PhoD-1101Si LP PhotoPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105Si C PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109Spheres fromSpheresSpheres fromTD-1102InGaAs PhoTD-1103Spheres from	otodiode	350-1100	-	5 Ø	30	W nm	11/30	91	82/83
LP-9920     InGaAs Pho       LP-9930     Si & InGaAs Pho       LPPA-9901     Si & InGaAs Pho       PD-1101     Si LP Photo       PD-1102     InGaAs Pho       PD-1103     InGaAs Pho       PD-1104     InGaAs Pho       PD-1105     SiC Photo       PD-1106     Si LP Photo       PD-1107     Si LP Photo       PD-1108     Si LP Photo       PD-1109     GaP Photo       PD-1109     Spheres from       Spheres     Si Photoc       TD-1101     Si Photoc       TD-1102     InGaAs Pho       Integrating     Spheres from       TD-1103     InGaAs Pho       TD-1103     InGaAs Pho	hotodiode	400-1100	-	2 Ø	8	W nm	11/8	91	82/83
LP-9930Si & InGaAs PLPPA-9901Si LP PhotoPD-1101Si LP PhotoPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105Si C PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109Spheres fromSpheresSpheres fromTD-1101Si PhotoTD-1103InGaAs PhoIntegratingSpheres fromSpheresSpheres from	otodiode	400-1100	-	12.5 Ø	50	W nm	11/50	91	83
LPPA-9901 Si LP Photo PD-1101 Si Photo PD-1102 InGaAs Pho PD-1103 InGaAs Pho PD-1103 Si C Photo PD-1105 Si C Photo PD-1106 Si LP Photo PD-1106 Si LP Photo PD-1107 Si LP Photo PD-1108 Si LP Photo PD-1109 GaAs Pho Integrating Spheres Since Strom Spheres Since Strom TD-1101 Si Photoco TD-1102 InGaAs Pho InGaAs Pho InGaAs Pho	Photodiode	800-1800	-	12.5 Ø	50	W nm	11/50	91	83
PD-1101Si PhotocPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotocTD-1102InGaAs PhoIntegrating Spheres from	s Photodiode	400-1800	-	12.5 Ø	50	W nm	11/50	91	83
PD-1101Si PhotocPD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotocTD-1102InGaAs PhoIntegrating Spheres from		Laser Pu	ulse Shape I	Radiant Powe	er Detector	<u> </u>			
PD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotoTD-1102InGaAs PhoIntegrating SpheresSpheres from	hotodiode	400-1060	-	12.5 Ø	50	V / W nm	LPPA	91	84
PD-1102InGaAs PhoPD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotoTD-1102InGaAs PhoIntegrating SpheresSpheres from		Co	mpact Size	Modular Dete	ectors		1		
PD-1103InGaAs PhoPD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhotoPD-1110GaAsP PhotoIntegrating SpheresSpheres fromTD-1101Si PhotocoTD-1102InGaAs PhoIntegrating SpheresSpheres from	otodiode	200-1100	-	3.6 x 3.6	-	With Accessory	11	90	85
PD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotoTD-1102InGaAs PhoIntegrating SpheresSpheres from	Photodiode	800-1800	-	0.3 Ø	_	With Accessory	11	90	85
PD-1104InGaAs PhoPD-1105SiC PhotoPD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1109GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotoTD-1102InGaAs PhoIntegrating SpheresSpheres from	Photodiode	800-1800	-	1Ø	_	With Accessory	11	90	85
PD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1110GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotocTD-1102InGaAs PhoIntegrating SpheresSpheres from	Photodiode	800-1800	_	3Ø	_	With Accessory	11	90	85
PD-1106Si LP PhotoPD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1110GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotocTD-1102InGaAs PhoIntegrating SpheresSpheres from	notodiode	200-400	_	1.1 x 1.1	_	With Accessory	11	90	85
PD-1107Si LP PhotoPD-1108Si LP PhotoPD-1109GaP PhotoPD-1110GaAsP PhoIntegrating SpheresSpheres fromTD-1101Si PhotocoTD-1102InGaAs PhoTD-1103InGaAs PhoIntegrating SpheresSpheres from		400-1050	_	0.5 x 0.5	_	With Accessory	11	90	85
PD-1108Si LP PhotoPD-1109GaP PhotoPD-1110GaAsP PhotoIntegrating SpheresSpheres fromTD-1101Si PhotocoTD-1102InGaAs PhotocoTD-1103InGaAs PhotocoIntegrating Spheres from		400-1050	_	1.9 x 1.9	_	With Accessory	11	90	85
PD-1109 GaP Photo PD-1110 GaAsP Photo Integrating Spheres from Spheres TD-1101 Si Photoco TD-1102 InGaAs Pho TD-1103 InGaAs Pho Integrating Spheres from		400-1050	-	3.4 x 3.4	_	With Accessory	11	90	85
PD-1110 GaAsP Pho Integrating Spheres Spheres from TD-1101 Si Photoc TD-1102 InGaAs Pho TD-1103 InGaAs Pho Integrating Spheres from		250-550	-	1.1 x 1.1	_	With Accessory	11	90	85
Integrating SpheresSpheres fromTD-1101Si PhotocoTD-1102InGaAs PhoTD-1103InGaAs PhoIntegratingSpheres from		200-760	_	2.3 x 2.3				+	
TD-1102InGaAs PhoTD-1103InGaAs PhoIntegratingSpheres from		1 inch) to 500 mr				an be Combined w	ith PD-11		ntegrating oheres
TD-1102InGaAs PhoTD-1103InGaAs PhoIntegratingSpheres from		Tempera	ature Contro	olled Modular	Detectors			· · · ·	
TD-1102InGaAs PhoTD-1103InGaAs PhoIntegratingSpheres from	otodiode	200-1100	-	3.6 x 3.6	-	With Accessory	TD11	90	86
TD-1103 InGaAs Pho Integrating Spheres from		800-1800	-	0.3 Ø	-	With Accessory	TD11	90	86
Integrating Spheres from		800-1800	-	1Ø	-	With Accessory	TD11	90	86
			ם (20 inch) ח		l fered which a	can be Combined w			ntegrating
•			or Radiant Po	wer Measurer					oheres
MD-37-GP6 GaP Photo	notodiode	250-550	-	2.5 x 2.5	-	With Accessory	MD37	90	86

# Page 62 - Light Detectors



# Table of Contents

Model	Spectral Function	Wavelength Range	Diffuser Diameter	Aperture (mm)	Sphere Diameter	Unit	Pack	age	Specifi- cation
		nm	mm	F.O.V.	mm		Туре	Page	Page
			Modula	r Detectors	4				
MD-37-GAP	GaAsP Photodiode	200-790	-	2.3 x 2.3	-	With Accessory	MD37	90	86
MD-37-IGA1	InGaAs Photodiode	800-1800	-	1Ø	-	With Accessory	MD37	90	86
MD-37-IGA3	InGaAs Photodiode	800-1800	-	3 Ø	-	With Accessory	MD37	90	86
MD-37-SC1	SiC Photodiode	200-400	-	1 x 1	-	With Accessory	MD37	32	86
MD-37-SU13	Si Photodiode	200-1100	-	3.6 x 3.6	-	With Accessory	MD37	32	86
MD-37-SU33	Si Photodiode	200-1100	-	5.8 x 5.8	-	With Accessory	MD37	32	86
MD-37-SU100	Si Photodiode	200-1100	-	10 x 10	-	With Accessory	MD37	32	86
MD-37Z-C11	Front Ca	ap Adapter		11 Ø	-	-	MD37	32	86
MD-37Z-DR11	RADIN Diffuser Cap	200-2500	-	11 Ø	-	-	MD37	32	86
MD-37Z-H7	Filter Adapter for 0.5 or 1	inch Diameter F	ilter				MD37	32	86
MD-37Z-H7xx	Filter Adapter with Custor	m Specified Ban	d-pass or We	edge Filter			MD37	32	86
Integrating Spheres	Spheres from 100 mm (4			ameter are Of wer Measurer		an be Combined wi	th MD-37		ntegrating oheres
		1 cm <sup>2</sup>	Circular Are	ea Modular D	etectors				
PD-9304	Si Photodiode	350-1100	-	11.28 Ø	-	W nm	PD93	91	87/88
PD-9306	UV-Si Photodiode	250-1100	-	11.28 Ø	-	W nm	PD93	91	87/88
PD-9302-VL	Photopic Front Filter	380-780	-	15.5 Ø	-	-	PD93	91	87/88
PD-9302-RW	Radiometric Front Filter	400-1100	-	11.28 Ø	-	-	PD93	91	87/88
PD-9303-COS	Diffuser Front Cap	380-1100	20	Cosine	-	-	PD93	91	87/88
		Avera	ged LED Int	ensity Front	Adapter			1	
PD-9311	0.01 Steradian Solid Ang	e Front Tube for	PD-9304; CI	E 127 Conditi	on B	cd or W/sr	PD93	91	87/88
PD-9312	0.001 Steradian Solid Ang	gle Front Tube fo	or PD-9304; (	CIE 127 Condi	tion A	cd or W/sr	PD93	91	87/88
			Weather	Proof Option					
WPD-01	Weather Proof Detector H	lousing to be Ec	quip with Pho	tometric or Ra	diometric Det	ector	WPD	92	88
			Light Detect	or Accessor	es				
ИМВ	Integrating Spheres from	100 mm / 4 in to	o 500 mm / 20	) in Diameter	with Accessor	ry for System Set-up	S	Int.	Spheres
UP	Precise Machined Integr	ating Spheres fr	om 8 mm / 0.	3 in to 10 in D	ia. with Acce	ssory for System se	t-ups	Int.	Spheres
SRT-M37	Front Lenses, Barrels & Adapters for 37-type Detectors 37 89							89	
/WQ	Waterproofing Option for 37-type Detectors 92							90	
/LC	Cable Length Extension								90
			Ca	rrying Cases					90

### **Optometer Model Depending Plug Types:**

Gigahertz-Optik light detectors come equipped with a plug type compatible with the particular optometer it will used with (see table).

The -2 calibration data plug contains an Eeprom that stores the detectors serial number, model, units and calibration data or data table for automatic information transfer to the meter when connected.

Both standard BNC (-1) and ITT-9PIN (-4) connectors are straight forward plugs without internal Eeprom.



- 1 BNC Connector

- 2 Calibration Data Connector

- 4 ITT Type Connector

Plug Type		Meter
BNC	-1	P-9202, P9801, P-9802, TR-9600
Calibration Data	-2	P-9710, P-2000
ITT (MSDM9)	-4	X9, X1, HCT-99

# Measurement Range Calculation; Illuminance & Irradiance

#### Calculation of the Maximum Measurable Intensity and Resolution:

Typical sensitivity factors for each light detector in its calibrated units are included in the specifications listed on the following pages. The maximum measurable intensity and attainable resolution with the different models of optometer offered can be calculated as follows

is

#### **Maximum Resolution**

= Meter Current Signal Resolution / Detector Sensitivity e.g. 0.1 pA / 3 nA / ( mW / cm<sup>2</sup> ) = 0.33 nW / cm<sup>2</sup> Maximum Signal = Max Meter Signal Current / Detector Sensitivity e.g.  $20 \,\mu\text{A} / 3 \,\text{nA} / (\,\text{mW} / \text{cm}^2\,) = 6.66 \,\text{W} / \text{cm}^2$ 

tor heads".

 $lx = lm / m^2$ 

 $fc = Im / ft^2$ 

Illuminance E

vertically in the units of lux (lx) or Illuminance refers to the luminous footcandles (fc). When the inciflux per unit area impinging on a particular surface from one or dent light is not parallel, which in more sources of light. Illuminance practice is normal, a cosine diffuser must be used as the measmeasured horizontally and

#### Irradiance E

Irradiance refers to the radiant power impinging on a particular surface from one or more sources of optical radiation. Irradiance is measured horizontally and verti-

Each optical radiation quantity

calls for its own measurement

geometry for the proper reception

of the radiation. Illuminance and

irradiance measurement systems require a cosine corrected field of cally in the units of watts per square meter (W/m<sup>2</sup>). When the incident optical radiation is not parallel, which in practice is normal, a cosine diffuser

**Cosine Corrected Field of View** 

view. Only then can the laws for

the incidence of diffuse radiation

from one or more sources be sat-

If the incident radiation is not par-

allel, the accuracy of the cosine

must be used to correct the measurement geometry for irradiance detector heads.

function is critically important to

the result of the measurement.

urement geometry for "lux detec-

VL-1101 with UPB-50-L

	$\Phi^*\cos \alpha$



isfied

Photopic Spectral Response

Ordering Inf	Ordering Information & typical Specifications										
Model	λresp	f <sub>1</sub>	Typical Sensitivity	Imax Sensing Area cable		cable	Temp.	plug	package		
	InGaAs	$\leq$ %	nA/lx mA 0.3 1		photodiode	photodiode m			page		
VL-1101	V(λ)     5     0.3     1     1.1 mm Ø     2     5 - 40     1,2,4										
K-VL1101-I	Calibratio	on of t	he illuminance sensitivi	ty in Ix	or fc. Including	K-SR in	order with	new de	etector		
KDW-S	Calibratio	on of t	he spectral sensitivity ir	n comb	pination with acc	essory					
KDW-P	Calibratio	on of t	he integral photometric	sensi	ivity in combinat	ion with	n accessory	/			

#### VL-1101: Compact Size Illuminance Detector



of discrete linear arrays (11 mm grid size) to measure the luminous intensity distribution of spot lamps.

The detector can be combined with integrating spheres with 11type adapter for luminous flux measurements.

Calibration is offered with or without accessory.

#### VL-3701, VL-3702, VL-3704 & VL-3705: Photopic & Scotopic Illuminance Detectors

Illuminance is one of the most important measurements in the field of commercial or industrial lighting. It enables us to see things properly and comfortably. Illuminance detector's spectral responsiveness must be equal the photopic vision characteristic  $V(\lambda)$ , as defined by the CIE. The very best matching of the detector spectral response to the photopic vision function is necessary especially if the measurements involve different kinds of light sources such as daylight, tungsten, fluorescent lamps, LED's etc

Illuminance detectors require a cosine field-of-view characteristic. This spatial response properly 'weights' the effects of the illumination reaching the detector from different angles of incidence. The accuracy of the cosine correction is as equally important as the photometric responsivity for illuminance detectors

The height of the detector is an other important factor, since the illuminance should be measured close to the subject surface. To measure the ANSI Lumen rating of projectors a max. height of 20 mm is specified.



#### VL-3701: Photopic with DIN Class A Parameters

The VL-3701 is a high quality detector for the measurement of illuminance. Its specifications satisfy the requirements for lux detector heads according to DIN 5032 Part 7.

The VL-3701's spectral adaptation to the  $V(\lambda)$  function for daylight vision is better than 3%. The precise adaptation of its spectral sensitivity to the ideal V( $\lambda$ ) function makes this detector head appropriate where sharply different spectra need to be measured (arc lamps, color monitors, LEDs, etc.). The deviation from the cosine function is < 1.5 %. The low detector profile of 20 mm allows its use in applications where the height of the detectors surface over the illuminated surface is restricted.

Calibration of the illuminance, in lux or fc is performed at Gigahertz-Optik's Calibration Laboratory for Optical Radiation Quantities, and is supported by a factory certificate.

#### VL-3702: Photopic with DIN Class B Parameters

The VL-3702 detector is constructed the same as the VL-3701 model, but has a lower precision of photometric adaptation, at 6 %, sufficient for DIN quality class B, and a maximum deviation of 3 % from the cosine function. Available versions and calibration are the same as model VL-3701.

#### VL-3704: Photopic with DIN Class B Parameters

The VL-3704 detector specifications are equal to that of the VL-3702 detector except that it is less sensitive. This allows the VL-3704 detector to be used with optometers like the X11, X91 and P-9802 with amplifiers having a limited signal range.

#### VL-3705: Scotopic

The spectral response of the human eye changes with light intensity. In daylight conditions >0.1 lx the photopic action spectrum applies. In low brightness applications 0.01 to 0.0001 lx the spectral eye response changes to the scotopic (dark adapted) action spectrum. Light intensities below 0.0001 k can not be detected by the human eye.

The VL-3705 detector offers a precise scotopic spectral match





Scotopic Spectral Response

with an  $f_1 = < 5\%$ . This is within DIN class B specifications qualifying the VL-3705 for use in most low light level application.

Ordering In	formation	& typ	ical Specif	ications						
Model	λresp	f <sub>1</sub>	Typical S	Sensitivity	Imax	Sensing Area	cable	Operation	plug	package
		$\leq$ %	nA/lx	nA/fc	mA	Diffuser	m	Temp.		page
VL-3701	$V(\lambda)$	3	0.5	5	1	7 mm Ø	2	5 - 40	1,2,4	90
VL-3702	$V(\lambda)$	5	0.5	5	1	7 mm Ø	2	5 - 40	1,2,4	90
VL-3704	$V(\lambda)$	5	20 pA/lx	200 pA/fc	1	7 mm Ø	2	5 - 40	4	90
VL-3705	V'(λ)	5	0.5	5	1	7 mm Ø	2	5 - 40	1,2,4	90
K-VL37xx-I	Calibratior xx = mode			ensitivity in ly	or fc.	Including K-SR	with ne	w detector	order.	
KDW-S	Calibratior	n of sp	ectral sensi	tivity with ac	cesso	y components				
KDW-P	Calibratior	n of int	egral photo	metric sens	itivity i	n combination w	ith acce	essory		
K-FOV	Calibratior	n of co	sine field o	f view functi	on. Inc	luded with new	detecto	or order		
V/L 07:04 III.					-000	Negel imite (8/	•			

VL-37xx Illuminance Detectors compa	arison DIN	I 5032 Class	Limits (%)			
Characteristic	Symbol	DIN 5032 Class A	DIN 5032 Class B	VL-3701	VL-3702	VL-3704
Calibration Uncertainty	$U_{kal}$	1.5	3	1.1	1.1	1.1
V(λ) Match Characteristic	f <sub>1</sub>	3	6	3	5	5
UV Response Characteristic	u	1	2	0.01	0.01	0.01
IR Response Characteristic	r	1	2	0.01	0.01	0.01
Directional Response (Cosine)	f <sub>2</sub>	1.5	3	1.5	3	3
Linearity Characteristic	f <sub>3</sub>	1	2	0.2	0.2	0.2
Fatigue Characteristic (at 1 klx)	f <sub>5</sub>	0.5	1	0.1	0.1	0.1
Temperature Dependence Characterist.	f <sub>6</sub>	2	10	1	1	1

#### **Accessories:**

Mating SRT-M37 front lens adapters allow VL-37 detectors to be adapted for luminance measurements. See page 89 for specifications.



VL-37 detectors can be combined with integrating spheres for luminous flux measurements. See integrating sphere chapter.





# Light Detectors - Page 65 Illuminance / Irradiance

### PD-9310: High Sensitivity Photopic Illuminance Detector



**Ordering Information & typical Specifications** 

 $f_1$ 

≤%

3

6

λresp

V(λ)

V(λ)

Typical Sensitivity

nA/lx

2.8

2.8

Calibration of the spectral sensitivity in combination with accessory

Imax

mΑ

1

1

Calibration of the integral photometric sensitivity in combination with accessory

Calibration of the photometric sensitivity in Ix or fc. Including K-SR in order with new detector

Sensing Area

Diffuser

30 mm Ø

30 mm Ø

cable

m

2

2

Operation

Temp

5-40 °C

5-40 °C

Model

PD-9310A

PD-9310B

K-PD9310-I

KDW-S

KDW-P



large 30 mm diameter diffuser window combined with a large 100 mm<sup>2</sup> area photodiode. It's high sensitivity makes this detector ideal for all applications in which very low light levels needs to be measured.

The detector can be combined with the PD-93 modular detector accessories and with integrating spheres with 37-type adapters. Calibrations are offered with or without accessory.



### TD-11VL01: Temperature Controlled Illuminance Detector

plug package

1,2,4

1,2,4

page

91

91



The TD-11VL01 detector offers a unique temperature control for filter and photodiode by precise heating to a higher operation temperature than the ambient temperature. The control electronic is supplied in an external housing with plug-in power supply. The detector can be combined with integrating spheres with 11-type adapter for luminous flux measurements.

Calibrations are offered with or without accessory.

Ordering Infor	Ordering Information & typical Specifications											
Model	λresp	f <sub>1</sub>	Typical Sensitivity	Imax	Sensing Area	cable	Opera-	plug	package			
		≤%	nA/lx	mA	photodiode	m	tion Temp.					
<b>TD-11VL01</b> V(λ)     5     tbc with sphere     1     2.77 mm x 2.77 mm     2     55°C     1,2,4     90												
K-TD11VL01-I	Calibra	tion of	the photometric sen	sitivity	in lx or fc. Including k	K-SR in	order with	new o	detector			
KDW-S	Calibra	tion of	the spectral sensitiv	ity in c	ombination with acce	ssory						
KDW-P	Calibration of the integral photometric sensitivity in combination with accessory											

VISIBLE-NIR Polychromatic Radiometric Light Detectors

Polychromatic radiometric detectors are primarily used for integral intensity measurements in a defined range of wavelengths. Correction filters are used to shape the bare detector re-



RW-01 BLUE 400-500 nm Spectral Response

sponse to the desired bandpass function, and ensure a radiometric responsiveness function. The plots below and on the following page summarizes the standard measurement heads

offered with correction filters and spectral responses. Customerspecific solutions with respect to spectral range, measurement geometry and other design fea-

tures are available along with full

assistance from our sales engineers.



sensitivity 0.5 ē 0 350 450 550 650 750 850 wavelength (nm)



# Irradiance

#### VISIBLE-NIR Polychromatic Radiometric Light Detectors



#### RW-11: Compact Size VISIBLE-NIR Polychromatic Radiometric Detectors



**The RW-11** detectors offer a compact package with flange mount which allows arrangement of discrete linear arrays with a grid of 11 mm to measure the radiant intensity distribution of light sources. The detector can be combined with integrating spheres with 11type adapter for radiant power measurements. Calibrations are offered with or without accessory.

RW-11 assembled to UM series sphere



🖲 Gigahertz-Optik

ordering in	ordering miorination & typical Specifications											
Model	λresp	Response	Typical Se	ensitivity	Imax	Sensing Area	cable	Operation	plug	package		
	polychromatic	Page	nA/W/m² nA/W		mA	photodiode	m	Temp.		page		
RW-1103	VISIBLE400-800 nm	65			1	1.1 x 1.1 mm	2	5-40 °C	1,2,4	90		
RW-1104	NIR <sub>800-1000 nm</sub>	66	tbc with	sphere	1	1.1 x 1.1 mm	2	5-40 °C	1,2,4	90		
RW-1105	VISNIR400-1000 nm	66			1	1.1 x 1.1 mm	2	5-40 °C	1,2,4	90		
K-RW11xx-I	Calibration of the irrad	iance sensiti	vity in A/W/m² a	and A/W/cm <sup>2</sup> .	Including	g K-SR in order wit	h new de	etector. xx = de	etector m	nodel		
KDW-S2	Calibration of the spec	Calibration of the spectral irradiance sensitivity at one or a chain of wavelength without or in combination with accessory										
KDW-R	Calibration of the integ	gral radiomet	ric sensitivity in	combination	with acce	essory						

# TD-11RW: Temperature Controlled VISIBLE-NIR Polychromatic Radiometric Light Detectors



**The TD-11RW** detectors offers a unique temperature control for filter and photodiode by precise heating to a higher operation temperature than the ambient temperature. The control electronic is supplied in an external housing with plugin power supply. The detector can be combined with integrating spheres with 11type adapter for radiant power measurements. Three different kinds of spectral responses e.g. TD-11RW03, TD-11RW04 and TD-11RW05 are offered which fulfill the need of most radiometric applications for integral detectors. Calibrations are offered with or without accessory.

Ordering Information	ation & typical Specif	ications									
Model	λresp	Response	51 5		Imax Sensing Area		cable	Operation	plug	package	
	polychromatic	Page			mA	photodiode	m	Temp.		page	
TD-11RW03	VISIBLE400-800 nm	65			1	2.8 x 2.8 mm	2	55°C	1,2,4	90	
TD-11RW04	NIR <sub>800-1000 nm</sub>	66	tbc with	sphere	1	2.8 x 2.8 mm	2	55°C	1,2,4	90	
TD-11RW05	VISNIR400-1000 nm	66			1	2.8 x 2.8 mm	2	55°C	1,2,4	90	
K-TD-11RWxx-I	Calibration of the irrad	iance sensiti	vity in A/W/m <sup>2</sup> a	and A/W/cm <sup>2</sup> .	Includ	ling K-SR in order	with nev	v detector			
KDW-S	Calibration of the spec	alibration of the spectral sensitivity in combination with accessory									
KDW-R	Calibration of the integ	gral radiomet	ric sensitivity in	combination	with a	ccessory					

Irradiance

# RW-37: High Sensitivity VISIBLE-NIR Broadband Irradiance Detectors



The RW-37 irradiance detectors are for use in radiometric applications requiring high sensitivity. The diffuser window offers a cosine corrected field-of-view. The compact 37-type package fits

into tight spaces in end-on viewing set-up. A side M6 tapped hole allows post mounting or fixturing into the application.

The detector can mounted to integrating spheres using 37type adapters and calibrated for radiant power measurements. See Integrating Spheres chapter. For simple radiometric radiance measurements the RW-37 detectors can be combined and calibrated with SRT-M37 front lenses. Calibrations are offered with or without components attached.



Ordering Info	rmation & typical Spe	cifications								
Model	λresp	Response	Typical Se	ensitivity	Imax	Sensing Area	cable	Temp.	plug	package
	Broadband	Page	nA/W/m <sup>2</sup>	mA/W/cm <sup>2</sup>	mA	Diffuser	m	°C		page
RW-3701	BLUE <sub>400-500 nm</sub>	65	70	0.7	1	11 mm Ø	2	5-40	1,2,4	90
RW-3702	RED <sub>700-800 nm</sub>	65	70	0.7	1	11 mm Ø	2	5-40	1,2,4	90
RW-3703	VISIBLE400-800 nm	65	1.7	0.017	1	11 mm Ø	2	5-40	1,2,4	90
RW-3704	NIR <sub>800-1000 nm</sub>	66	200	2	1	11 mm Ø	2	5-40	1,2,4	90
RW-3705	VISNIR400-1000 nm	66	400	4	1	11 mm Ø	2	5-40	1,2,4	90
K-RW37xx-I	Calibration of radiome	tric sensitivity	in A/W/m² an	d A/W/cm <sup>2</sup> . In	cludes K	-SR in new detect	or order	. xx = detect	or model	
KDW-S	Calibration of spectral	sensitivity at	one or multiple	wavelengths	without	or in combination	with acc	cessory com	oonents	
KDW-R	Calibration of integral I	adiometric s	ensitivity in con	nbination with	accesso	ory components.				

ULTRAVIOLET broadband radiometric detectors are primarily used for integral intensity measurements in defined UV band pass ranges of wavelengths.



lowing page summarize the standard measurement heads offered with correction filters and spectral responses. Customerspecific solutions with respect to

sensitivity 60

<u>e</u>

**ULTRAVIOLET** Polychromatic Irradiance Light Detectors

spectral range, measurement geometry and other design features are available along with full assistance from our sales engineers.



UV-01 Spectral Response: UV-A 315-400 nm



UV-10 Spectral Response: UV-A 320-400 nm



sensitivity 2'0 Ð C 260 275 290 305 320 elength ( nm )

UV-11 Spectral Response: UV-B 280-320 nm



UV-03 Spectral Response: UV-C 200/250-280 nm



UV-16 Spectral Response: UV-A 305-400 nm

# Irradiance

#### **ULTRAVIOLET Broadband Irradiance Detectors**



#### UV-11: Compact Size UV Broadband Irradiance Detectors



**UV-A/UV-B** are the most frequently used wavelength ranges in the qualification procedure of light sources used in photobiological and industrial application. Besides standard narrow band detectors, Gigahertz-Optik offers real broadband detectors. **The UV-11** detectors offer a compact package with flange mount to arrange detector arrays with only a grid size of 11 mm. Or they can be combined with integrating spheres using 11type adapters for radiant power measurements. Calibrations are offered with or without accessory components.

Ordering Info	ormation & typical S	pecification	S									
Model	λresp	Response	Typical Se	ensitivity	Imax	Sensing Area	cable	Temp.	plug	package		
	Broadband	Page	nA/W/m <sup>2</sup>	nA/W	mA	photodiode	m	°C		page		
UV-1101	UV-A 315-400 nm	67	tbc	tbc with	0.2	2.5 x 2.5 mm	2	5-40	1,2 or 4	90		
UV-1102	UV-B <sub>280-315 nm</sub>											
K-UV11xx-I	Calibration of radion	netric sensitiv	vity in A/W/m <sup>2</sup> a	and A/W/cm <sup>2</sup> .	Includes	K-SR in with new	detector	order. xx =	detector ma	odel		
KDW-S	Calibration of spectr	ibration of spectral sensitivity at one or a chain of wavelength without or in combination with accessory components.										
KDW-R	Calibration of integra	bration of spectral sensitivity at one or a chain of wavelength without or in combination with accessory components.										

#### **TD-11UV: Temperature Controlled UV Broadband Radiometric Detectors**

**The TD-11UV** detectors offer a unique temperature control feature where the filter and photodiode are heated to a higher operating temperature than the ambient temperature.

The control electronics are supplied in an external housing with a plug-in power supply. The detector can be combined with integrating spheres with 11-type adapter for radiant power measurements. Two different kinds of spectral responses, UV-A<sub>315-400 nm</sub> and UV-B<sub>280-315 nm</sub> are offered which fulfill the need of most radiometric applications for integral detec-

tors. Calibrations are offered with or without accessory components.



Ordering Inform	nation & typical Spe	cifications										
Model	λresp	Response	Typical Se	ensitivity	Imax	Sensing Area	cable	Operation	plug	package		
	Broadband	Page	nA/W/m <sup>2</sup>	nA/W	mA	photodiode	m	Temp.		page		
TD-11UV01	UV-A 315-400 nm	67	tbc	tbc with	0.2	2.5 x 2.5 mm	2	55°C	1,2,4	90		
TD-11UV02	UV-B <sub>280-315 nm</sub>	67	tbc	sphere	0.05	1.1 x 1.1 mm	2	55°C	1,2,4	90		
K-TD11UVxx-I	Calibration of radiom	etric sensitivi	ty in A/W/m² ar	nd A/W/cm². I	ncludes	K-SR in new deteo	ctor ord	ler. xx = dete	ector			
KDW-S	Calibration of spectra	ibration of spectral sensitivity at one or multiple wavelengths without or in combination with accessory components										
KDW-R	Calibration of integra	Calibration of integral radiometric sensitivity with accessory components										

UV-37

### UV-37: High Sensitivity Broadband UV Irradiance Detectors

#### UV-3701, UV-3702, UV-3703

These detectors represents the classic DIN and CIE ultraviolet radiation wavelength bands UV-A\_{315-400 nm}, UV-B\_{280-315 nm} and UV-C\_{100-280 nm}. Deep UV in the range of 100 to 250 nm is very rarely

observed in radiometry. So the UV-3703 offers a practical UV-C solution with a sensitivity range from 250 to 280 nm.

#### UV-3710, UV-3711

Some non-biological disciplines specify the wavelength ranges 320 to 400 and 280 to 320 nm for the evaluation of UV-A and UV-B

eye operation).

tion >400 nm..

#### UV-3716

The sharply increasing shortwave spectral component of modern lamps creates a potential risk of exposure of the eyes to excessive UV-A radiation. The limit values for  $E_{\rm 305-400\ nm}$  irradi

In applications where the UV-A irradiance must measured amidst intense blue & visible light, the UV-3701 and UV-3710 may not offer sufficient visible wavelength range signal rejection. This will produce an inaccu-

UV-3717 JV-A rately high UV-A reading. A good example of this is fluorescent lamps which generate significant mercury lines at 405, 436, 546 and 577 nm. The UV-3717 is visi-

ble blind for best stray light rejec-

ance were established to limit

the intensity this radiation on the

retina as when optical fibers are

used in eye operations (aphakic



#### UV-3718

The model UV-3718 is designed for narrow band UV-light sources with distinct 254 nm peak emission spectrum. The RADIN type diffuser offers high UV stability under UV-C irradiation. The detector is calibrated with a narrow band light source at 254 nm.

#### UV-3719

The model UV-3719 offers a broadband spectral sensitivity from 250 to 390 nm. This detector is typically used in applications where the lamp emission spectrum is known allowing an a (z) calibration correction.

The UV-37 detectors are high sensitivity irradiance detectors for radiometric applications involging polychromatic optical radiation. The diffuser window provides a cosine corrected field-of-view for diffuse light signals. The compact 37-type package fits into tight spaces with a vertical end-on view. An M6 tapped hole allows post mounting or fixturing in the application. The detectors can be combined and calibrated with integrating spheres or SRT-M37 quartz front lenses.

Ordering In	formation & typical S	pecification	8							
Model	λresp	Response	Typical Se	nsitivity	Imax	Sensing Area	cable	Temp.	plug	package
	Broadband	Page	nA/W/m <sup>2</sup>	µA/W/cm <sup>2</sup>	mA	Diffuser	m	°C		page
UV-3701	UV-A 315-400 nm	67	4	40	0.1	11 mm Ø	2	5-40	1,2 or 4	90
UV-3702	UV-B <sub>280-315 nm</sub>	67	1.7	17	0.05	11 mm Ø	2	5-40	1,2 or 4	90
UV-3703	UV-C 200/250-280 nm	67	2	20	0.05	11 mm Ø	2	5-40	1,2 or 4	90
UV-3710	UV-A 320-400 nm	67	7	70	0.1	11 mm Ø	2	5-40	1,2 or 4	90
UV-3711	UV-B 280-320 nm	67	3	30	0.05	11 mm Ø	2	5-40	1,2 or 4	90
UV-3716	UV-A 305-400 nm	67	10	100	0.1	11 mm Ø	2	5-40	1,2 or 4	90
UV-3717	UV-A 325-400 nm	68	1	10	0.1	11 mm Ø	2	5-40	1,2 or 4	90
UV-3718	UV-C 254 nm	68	2	20	0.05	11 mm Ø	2	5-40	1,2 or 4	90
UV-3719	UV-A 250-400 nm	68	3	30	0.1	11 mm Ø	2	5-40	1,2 or 4	90
K-UV37xx-I	Calibration of radiome	etric sensitivity	in A/W/m <sup>2</sup> and	I A/W/cm². In	cludes k	<-SR in new detect	or order.	xx = detecto	or model	
KDW-S2	Calibration of spectra	l irradiance se	ensitivity at one o	or multiple wa	welengt	hs without or in co	mbinatio	n with acces	sory comp	onents
KDW-R	Calibration of radiome	etric sensitivity	with accessory	components	;					

#### UV-37: High Sensitivity Photobiologically Actinic UV Irradiance Detectors

Photobiology is concerned with the interaction of optical radiation and living organisms. Photobiological radiometric measurement systems are essential tools used for the evaluation of sources of optical radiation in cosmetics, medical and therapeutic, health and safety and other areas. The photodetectors used in these various application areas require spectral sensitivities corresponding directly to a particular photobiological response function. This enables direct measurement of the photobiological dose of a radiation source.

The spectral response function of the sensor is generated by optical correction filters used to tailor the bare detector response. Since the ideal target function cannot be perfectly duplicated, it is necessary to find the best possible compromise through extensive computer simulation taking a great variety of lamp spectra into account to minimize error. The associated series of measurements, which are an essential part of this engineering process, are carried out at Gigahertz-Optik's Calibration Laboratory for Optical Radiation Quantities. On the following page, irradiance light detectors available from Gigahertz-Optik that are primarily intended for applications in the fields of radiation protection and of radiation therapy are presented. For a specific response function not shown, or if you should require more detailed information, please contact our technical sales department directly.

### Actinic Irradiance

#### UV-3704: UV-Erythemal Actinic Irradiance Detector

The typical symptom of UV erythema is acute skin inflammation caused by UV radiation (sunburn) which has been linked to skin cancer over prolonged exposure times.

The UV-3704 detector permits direct measurement of the UV Erythema effective irradiance by means of it's spectral sensitivity adaptation to the nominal CIE Erythemal spectral function.

#### UV-3706: Bilirubin Actinic Irradiance Detector

Newborn jaundice or neonatal hyperbilirubinaemia, a yellowish appearance of the skin and whites of the eyes, is present to some degree in almost all newborn infants. This is caused by an elevated level of bilirubin molecule in the blood which results from immaturity of the liver function combined with the destruction of red blood cells present. When these levels are very high, one method of clearing the jaundice is by exposing the newborn to light in the blue spectral region between 400 to 550 nm. The UV-3706 detector allows the direct measurement of phototherapeutically effective irradiance employed for lowering the bilirubin level in newborn infants.

The UV-3708 detector permits direct measurement of effective

irradiance according to ACGIH/

Measurement uncertainty can be

evaluated and reduced using

Gigahertz-Optik's a(z) method

when the lamp emission spectra

is known with dataset provided.

ICNIRP guidelines.

Measurement uncertainty can be

evaluated and reduced using

Gigahertz-Optik's a(z) method

when the lamp emission spectra

is known with dataset provided.

#### UV-3708: ACGIH/ICNIRP Actinic Irradiance Detector

range.

The spectral weighting function for the acutely harmful effects of UV radiation, was developed by the American Conference of Governmental Industrial Hygienists (ACGIH) and the International Commission on Non-Ionising Radiation Protection (ICNIRP).

The spectral effectiveness in the UV-C and UV-B ranges is very high as compared to the UV-A

# UV-3709: Blue-Light Hazard Actinic Irradiance Detector

If optical radiation in the wavelength range between 380 and 550 nm of sufficient intensity reaches the retina it can cause photochemical injury. Radiation in the "blue" part of the spectrum triggers photochemical reactions, if the photon energy in the radiation is high enough, converting chemically unstable molecules into one or more other molecule types. BLH irradiance

is relevant exclusively in connection with ophthalmic instruments (microscopes, observation devices used in optician's practices). Otherwise the measurement of radiance applies. The UV-3709 detector head is designed for direct measurement of the effective irradiance for blue light hazard, responsible for retinal damage.









Ordering Int	formation & Specific	ations										
Model	λresp	Wavelength	Typical Se	ensitivity	Imax	Sensing Area	cable	Temp.	plug	package		
	Biological Effective	Range	nA/W/m <sup>2</sup>	$\mu$ A/W/cm <sup>2</sup>	mΑ	Diffuser	m	°C		page		
UV-3704	Erythema	250-400 nm	3	30	0.05	11 mm Ø	2	5-40	1, 2, 4	90		
UV-3706	Billirubin	400-550 nm	200	2000	1	11 mm Ø	2	5-40	1, 2, 4	90		
UV-3708	ACGIH/ICNIRP	250-400 nm	1,5	15	0.05	11 mm Ø	2	5-40	1, 2, 4	90		
UV-3709	Blue-light Hazard	400-520 nm	150	15000	1	11 mm Ø	2	5-40	1, 2, 4	90		
K-UV37xx-I	Calibration of actinic	irradiance sensitivit	y in A/W/m² an	ıd A/W/cm². Ir	ncludes	s K-SR in new dete	ector or	der. xx = det	ector			
KDW-S	Calibration of spectra	ation of spectral sensitivity with accessory components										
KDW-R	Calibration of radiom	ation of radiometric sensitivity with accessory components										



# Plant Physiology

#### PS-37: Plant Physiology Actinic Irradiance Detectors

The radiation conditions used in determining the rate of photosynthesis and the photosynthetic potential of various plant or algae types are not the same in all research institutions. Results obtained under very different radiation conditions, using detector heads with non-uniform rectangular (radiometric) characteristics, and then relating them to one another, may lead to false conclusions. This is because the varying spectra of the radiation sources in use are ignored in obtaining the measurement.

The solution is to evaluate the irradiance with a sensor with an appropriate spectral response function. It is presently assumed that the number of light quanta absorbed is responsible for plant growth, which implies that it is quantum magnitudes effective in plant biology that need to be measured. The most important

magnitude is the photosynthetic photon irradiance E  $_{\rm p,sy.}$ 

The wavelength limits of integration band pass need to be specified. If, for example, the photosynthetic photon irradiance is to be measured in the range of wavelengths between 400 nm and 700 nm, 320 nm and 500 nm and 590 nm to 900 nm, the integration is carried out in the corresponding spectral segments.

sensitivity

ē

0.5

0

380

430

480

530

580

wavelength (nm)

630

680

This numerical integration can be performed implicitly by means of a cosine corrected integral detector such as the PS-37 series.



Waterproof version of the PS-37 type detector head with dome cover

#### **PS-3701: PAR Irradiance Detector**

Photosynthesis is one of the most important biochemical processes on the planet. In the process of photosynthesis green plants absorb carbon dioxide from the atmosphere and water from the soil, combining them with the aid of radiation energy to build sugar, releasing oxygen and water into the atmosphere. In general plant physiology, the term **P**hotosynthetically **Active**  **R**adiation (**PAR**) refers to the radiation in the range of wavelengths between 400 nm and 720 nm.

The PS-3701 detector is designed for direct evaluation of **P**hotosynthetically **Active R**adiation in the 400 to 700 nm wavelength range. This detector can also be supplied with an optional water proof dome cover.

#### PS-3702: BLUE<sub>320-500 nm</sub> Plant Physiology Irradiance Detector

**Phototropism** describes the effect of optical radiation on the direction of plant growth. The regions of maximum effect lie in the blue range between 380 nm and 520 nm. Radiation can have the effect of causing parts of plants to move.

The PS-3702 detector is designed for direct evaluation of the **phototropism** effective irradiation over the spectral range from 320 to 500 nm for UV-A to green/visible assessment. This detector can also be supplied with an optional water proof dome cover.



#### PS-3703: RED<sub>560-900 nm</sub> Plant Physiology Irradiance Detector

Photomorphogenesis describes the way in which plants are formed under the influence of optical radiation. Radiation in the red region of the spectrum encourages linear growth, while blue radiation yields small, strong plants. To be more precise, the ratio of the radiation intensities in the range of wavelengths from 690 nm to 780 nm (long wavelength red) to the range of wavelengths from 560 nm to 680 nm (short wavelength red) is of great importance for the plant's biological processes. The PS-3703 detector evaluates photomorphogenesis effective radiation in the wavelength range from 590 to 900 nm. This detector can also be supplied with an optional water proof dome cover.



Ordering In	formation & typical \$	Specifications							
Model	λresp	Wavelength	Typical Sensitivity	Imax	Sensing Area	cable	Operation	plug	package
	Plant Physiology	Range	nA / µMol*m <sup>-2</sup> *s <sup>-1</sup>	mA	Diffuser	m	Temp.		page
PS-3701	PAR	400-700 nm	100	1	11 mm Ø	2	5-40°C	1, 2, 4	90
PS-3702	BLUE	320-500 nm	2	1	11 mm Ø	2	5-40°C	1, 2, 4	90
PS-3703	RED	590-900 nm	100	1	11 mm Ø	2	5-40°C	1, 2, 4	90
K-PS37xx-I	Calibration of PAR s	ensitivity in A/µMol*	m <sup>-2</sup> *s <sup>-1</sup> . Includes K-SR in n	ew deteo	ctor order. xx = de	etector n	nodel		
KDW-S	Calibration of spectra	al sensitivity at one o	or multiple wavelengths wit	hout or i	n combination wit	n access	sory compon	ents	
KDW-R	Calibration of integra	l radiometric sensiti	vity in combination with ac	cessory	components				

# High Intensity Irradiance / UV-Curing

**Three Different Package Designs** 

#### High Intensity Irradiance Detectors for UV-Curing Applications

Gigahertz-Optik offers irradiance detectors specially designed for hostile ambient conditions involving high intensity irradiation and high temperature.

The detectors consist of two

main components, the passive

RADIN element and the detector

The RCH type detectors are of-

that are connected by a flexible or rigid light guide.

The light guide protects the detector and corrective band pass filter from heat damage and also reduces measurement errors due to the temperature coefficient (drift) of the photodiode.

for use in the most common UV

RADIN itself is high UV irradiation and temperature stable up to 100°C with short peak measurements to 200°C. The low profile (9 mm) RADIN sensor element permits irradiance measurement close to the sample surface of the probe and offers a cosine

adapted field-of-view. Irradiances of up to 40 W/cm<sup>2</sup> can be measured.

For spot curing applications, adapters are available for simple positioning of different size light guide nozzles in front of the RADIN sensor.

#### **Three Different Spectral Sensitivities**

Detector spectral responses are

365 nm peak UVA as well as UV broadband and BLUE ranges.

> 410 430

450 470

i san il

06: UV Broadband Sensitivity







High Intensity & UV-Curing Actinic Irradiance / 360 Deg. F.O.V. Irradiance

### RCH-1xx, RCH-1xx & RCH-5xx: High Intensity Irradiance Detectors

Ordering Inf	ormation & typic	cal Spe	cifications	\$									
Model	λresp	Resp.	Max. E	E & T	Typ. Sensitivity	Imax	Sensing Area	Light Guide	cable	Front ad	dapter	plug	package
	Broadband	Page	W/cm <sup>2</sup>	°C	nA/mW/cm <sup>2</sup>	mA	Diffuser	cm / inch	m	Dia.	Page		page
RCH-006	UV Broadband	14	40	100	0.3	0.1	9 mm Ø	50 / 20	0.5	5,7,10	91	1,2,4	91
RCH-008	UVA Peak 365 nm	14	40	100	0.3	0.1	9 mm Ø	50 / 20	0.5	5,7,10	91	1,2,4	91
RCH-009	BLUE Peak 430 nm	14	40	100	0.2	0.1	9 mm Ø	50 / 20	0.5	5,7,10	91	1,2,4	91
RCH-106	UV Broadband	14	40	100	0.3	0.1	9 mm Ø	22 / 8.7	0.5	5,7,10	91	1,2,4	91
RCH-108	UVA Peak 365 nm	14	40	100	0.3	0.1	9 mm Ø	22 / 8.7	0.5	5,7,10	91	1,2,4	91
RCH-109	BLUE Peak 430 nm	14	40	100	0.2	0.1	9 mm Ø	22 / 8.7	0.5	5,7,10	91	1,2,4	91
RCH-506	UV Broadband	14	40	100	0.065	0.05	6 mm Ø	-	2	-	-	1,2,4	91
RCH-508	UVA Peak 365 nm	14	40	100	0.096	0.1	6 mm Ø	-	2	-	-	1,2,4	91
RCH-509	BLUE Peak 430 nm	14	40	100	0.065	0.1	6 mm Ø	-	2	-	-	1,2,4	91
K-RCHxxx-I	Calibration of irra	diance	sensitivity i	n A/W/a	m². Includes K-S	R in n	ew detector ord	er. xxx = dete	ector m	nodel			
KDW-S2	Calibration of spe	ectral irra	adiance se	nsitivity	at one or multiple	e wave	elengths without	or in combin	ation w	vith acces	sory c	ompor	nents
KDW-R2	Calibration of inte	egral irra	idiance ser	nsitivity	with accessory co	ompor	ients.						

#### ROD-360-UV18: 360 Degree F.O.V. UV-C<sub>254nm</sub> Irradiance Detectors

The ROD-360 is an unique irradiance detector which offers a 360 degree field-of-view. A quartz-rod is used to collect all

light irradiating it's diffuse detec-

(



tion window independent of the horizontal incident angle within the round angle. The vertical axis exhibits a diffuse viewing characteristic.

> With it's narrow band UV-C response the ROD-360 is suitable for measurement of the effective UVGI in air and water germicidal applications employing low and medium pressure mercury lamps.

Calibration is done using a 254 nm low pressure mercury light source.

The ROD-360 features a waterproof housing which allows measurements in humid or underwater applications. The 10 mm diameter clear quartz tube not only seals the detector rod but along with it's stainless steel housing the probe and active

ROD-360 fixed in vertical mound adapter with tripod threaded hole at the bottom

window can be easily cleaned. This makes the ROD-360 usable in gray water or other dirty measurement environments and also in medical applications.

A protective cap and mounting adapter is supplied to fixture the probe for vertical or horizontal use on a standard tripod or for integration into the application.

Ordering Information	on & typical Spe	ecifications									
Model	λresp	Wavelength	Typical S	ensitivity	Imax	Sensing Area	cable	Operation	plug	package	
	Narrow Band	Range	nA/W/m <sup>2</sup>	µA/W/cm <sup>2</sup>	mA	Rod	m	Temp.		page	
ROD-360-UV18	UV-C	220-280 nm	0.8	8	0.1	5 Ø x 6.5 mm	2	0-40°C	1,2,4	92	
K-ROD360-UV18-S	Calibration of th	ibration of the irradiance sensitivity in A/W/m <sup>2</sup> and. A/W/cm <sup>2</sup> Including K-SR in order with new detector.									
KDW-S2     Calibration of spectral irradiance sensitivity at one or multiple wavelengths											

# Multi-Channel Illuminance & Irradiance

### XD-95: Single & Dual Channel Illuminance & Irradiance Detector

The XD-95 detectors are designed for use with Gigahertz-Optik's multi-channel X1 and X9 hand-held optometers. The ability to measure two wavelength ranges without having to change and replace the detector head greatly simplifies the measurement procedure. Several detectors are offered covering common applications. A short overview describing the various models is shown here. More information is available in the X9 and X11 optometer data descriptions.

The XD-95 detectors are supplied in a compact lightweight

#### XD-9501: UV-A & UV-B Irradiance

Ultraviolet radiation is used in the treatment of certain skin diseases like Psoriasis and for photobiological studies like SPF testing. The phototherapist needs to monitor the source irradiance or energy for accurate dosimetry. Typical wavelength ranges of interest are UV-A<sub>320</sub>.

400nm, UV-B<sub>280-320nm</sub> and UV-B<sub>311nm</sub>. The XD-9501 detector offers a real polychromatic radiometric spectral match to the UV-A and UV-B spectral ranges and a precise cosine corrected field of view.

#### XD-9502: UV-A Irradiance & V(λ) Illuminance

Liquid Penetrant Testing of surface defects using the dye penetration process (DIN EN 1956, ASTM and MIL Standards) and Photostability testing of drug and drug products as stated in ICH guidelines (International Conference for Harmonization) requires the intensity or dose control of the light sources used for test sample irradiation. The XD-9506 detector offers a real broadband radiometric UV-A\_{315-400nm} detector and a illuminance detector (f1'  $\leq$ 6%), both with cosine diffuser for precise readings in both applications.

#### XD-9503: UV-B<sub>240-340nm</sub> & UV-A<sub>315-400nm</sub> Irradiance

TL01 and TL12 light sources are widely used in phototherapy applications. The phototherapeutic effective radiation of the TL01 and TL12 source falls within the 240 to 340 nm bandpass with peak at 313 nm.

tic effective UV-B and UV-A irradiance of TL01 and TL12 light sources. The detectors spectral sensitivity functions and calibration using TL1 and TL12 light sources ensure precise readings.

The XD-9503 is specially designed to measure the therapeu-

#### XD-9504: VISIBLE<sub>400-800nm</sub> & NIR<sub>800-1200nm</sub> Irradiance

Light has therapeutic effects in many different wavelength ranges. However broadband light sources will emit light at wavelengths outside the phototherapeutic wavelength range. This 'stray light' can interfere with the treatment or even cause harm in terms of heat for example, so must be measured in order to be kept to minimal levels.

The XD-9504 is designed to measure irradiance in the 'signal' wavelength range (400-800 nm) and 'noise' range (800-1200 nm). Both detectors offer a 'flat' broadband spectral response

#### XD-9505: UV-Broadband<sub>240-400nm</sub> Irradiance

In some applications the total UV irradiance of broadband light sources needs to be measured. But because of the unavailability of ultraviolet detectors with radiometrically 'flat' spectral sensitivities across the entire UV range, relative measurements are normal. The XD-9505 can be

used for relative broadband UV assessment and can provide more precise UV measurement through calibration correction when the emission spectrum of the light source is known.

The single detector offers a cosine diffuser for precise readings in diffuse irradiation zones. and low profile plastic housing with a 2 m long cable with -4

type plug for use with the X9 and X1 meters.





wavelength ( nm )











# 🕞 Gigahertz-Optik

# Light Detectors - Page 75





# **Multi-Channel Illuminance & Irradiance**

#### XD-9506: EM6 & EN14255 ACGIH Irradiance

Several guidelines and standards concerning UV hazard risk assessment exist today. Most often the ACGIH/ICNIRP spectral weighting function for the acutely harmful effects is recommended as reference for irradiance measurements (American Conference of Governmental Industrial Hygienists (ACGIH) and the International Commission on Non-Ionising Radiation Protection (ICNIRP). The XD-9506 ACGIH detector is designed to meet EM 6, DIN-EN 14255 and BGV B9 regulations. It's two cell design offers best possible isolation between UV-A and UV-B.

As part of these requirements an additional UV-A irradiance detector like the XD-9502, XD-9511 or UV-3717 is needed.

#### XD-9509: 2-Channel Erythema Effective Irradiance

The typical symptom of UV erythema is acute skin inflammation caused by UV radiation (sunburn) which may be linked to skin cancer over prolonged exposure times. The XD-9509 detector permits direct measurement of the UV Erythemal effective irradiance by means of it's spectral sensitivity adaptation to the nominal CIE Erythemal spectral function. It's two cell design offers best possible isolation between the UV-A and UV-B range.

#### XD-9510: DIN EN 12198 ACGIH Irradiance

The XD-9510 detector is a modified version of the XD-9506 ACGIH irradiance detector. The difference is the spectral response of the UV-A cell which simulates the *flat UV-A response* specification for DIN EN 12198. Both detectors offer a cosine diffuser.

# XD-9511: UV-A & EN 14255 UV-A<sub>315-400nm</sub> Irradiance

The XD-9511 offers a single detector with UV-A spectral response from 315 to 400 nm. Because of it's good blocking of the visible spectrum above 400 nm, this detector can also be used in combination with the XD-9506 for measurement of UV-A irradiance following DIN EN

#### 14255 and BGV B9.

The detector offers a cosine diffuser for precise measurements in diffuse light incident conditions.

Ordering In	formation & typ	oical Specificati	ons								
Model	Ch1 λresp	Ch2 λresp	Typical Se	ensitivity	Sensing Area	cable	Operation	Optor	neter	plug	package
			nA/	nA/	Diffuser	m	Temp.	X9 type			page
XD-9501	UVA <sub>320-400nm</sub>	UVB <sub>280-320nm</sub>	17 W/m <sup>2</sup>	2.5 W/m <sup>2</sup>	8 mm Ø	2	5-40°C	X96	X11	4	91
XD-9502	UVA <sub>315-400nm</sub>	$V(\lambda)$	10 W/m <sup>2</sup>	0.05 lx	8 mm Ø	2	5-40°C	X98	X11	4	91
XD-9503	UVA 315-400nm	UV <sub>Peak 311nm</sub>	4 W/m <sup>2</sup>	2 W/m <sup>2</sup>	8 mm Ø	2	5-40°C	-	X11	4	91
XD-9504	VIS <sub>400-800nm</sub>	NIR <sub>800-1200nm</sub>	0.7 W/m <sup>2</sup>	130 W/m <sup>2</sup>	8 mm Ø	2	5-40°C	-	X11	4	91
XD-9505	UV <sub>Broadband</sub>	-	13 W/m <sup>2</sup>	-	8 mm Ø	2	5-40°C	X97	X11	4	91
XD-9506	UVA <sub>ACGIH</sub>	UVB <sub>ACGIH</sub>	4000 W/m <sup>2</sup>	1.5 W/m <sup>2</sup>	8 mm Ø	2	5-40°C	-	X11	4	91
XD-9509	UVA <sub>Erythema</sub>	UVB <sub>Erythema</sub>	tbc W/m <sup>2</sup>	tbc W/m <sup>2</sup>	8 mm Ø	2	5-40°C	-	X11	4	91
XD-9510	UVA <sub>ACGIH</sub>	UVB <sub>ACGIH</sub>	4000 W/m <sup>2</sup>	1.5 W/m <sup>2</sup>	8 mm Ø	2	5-40°C	-	X11	4	91
XD-9511	UVA <sub>315-400nm</sub>	-	4 W/m <sup>2</sup>	-	8 mm Ø	2	5-40°C	X97	X11	4	91
K-XD95xx-l	Calibration of int	tegral illuminance	e/irradiance sen	sitivity in A/Ix	or A/W/cm <sup>2</sup> . Inc	ludes ł	<-SR in new	detector o	rder. xx =	detecto	r model
KDW-S2	Calibration of sp	oectral irradiance	sensitivity at or	e or multiple	wavelengths						

# Multi-Channel Illuminance, Irradiance & Luminous Color

### **TP-4501: Triple PAR Irradiance Detectors**



The TP-4501 is a unique detector head designed for use with the X11 hand-held optometer to measure on-site lighting conditions in greenhouses, plant growth studies, soil and seed science.

The compact detector houses

four different detectors:

PAR: Measures the photonsynthetically effective radiation in the 400 to 700 nm wavelength range actinic for green growth. BLUE: Measures the phototropism effective radiation within the 320 to 500 nm wavelength range actinic for directional plant growth.

RED: Measures the photomorphogensis effective radiation within the 560 to 900 nm wavelength range actinic for plant growth formation.  $V(\lambda)$ : Measures the **photomet-**

ric illuminance which specifies the ambient lighting condition. The low profile detector housing



is only 20 mm high with a central 10 mm diameter diffuser which ensures the cosine corrected field of view required in illuminance & irradiance applications.

Ordering In	formation &	Specification	ns											
Model	λresp	λresp	λresp	λresp	Т	ypical S	Sensitivi	ty	0	cable	Operation			package
	Ch1	Ch2	Ch3	Ch4	CH1	CH2	CH3	CH4	Area	m	Temp. ℃	Туре		page
	photometric	Plant I	Physiology Eff	ective	pA/lx	pA/μ	Mol*m	<sup>-2*</sup> S <sup>-2</sup>	Diffuser		C			
TP-4501	$V(\lambda)$	PAR <sub>400-700 nm</sub>	BLUE <sub>320-500 nm</sub>	RED <sub>560-900 nm</sub>	40	0.1	2	70	10 mm Ø	2	5-45	X11	4	90
K-TP4501-I	Calibration o	f photometric	& plant physic	logy effective	e sensit	ivity in	A/Ix & A	õMol*	m <sup>-2*</sup> s <sup>-2</sup> . Inc	ludes l	≺-SR in new	detect	or ord	er.
KDW-S2	Calibration o	f spectral irrad	diance sensitiv	ity at one or r	nultiple	wavele	ngths v	vithout (	or in combi	nation	with access	sory cor	npone	ents.

#### CT-3701 & CT-4501: Luminous Color and Illuminance Detectors



#### CT-3701

The CT-3701 and CT-4501 luminous color detectors are designed with precise x,y,z tristimulus matching functions to measure the illuminance, x,y chromaticity values and color temperature of broad band emitting light sources. Four detectors mounted in one compact housing are used to form the  $x_{\mbox{\tiny short}},\,x_{\mbox{\tiny long}},\,y$  and z spectral tristimulus values for the 2 degree viewer.

CT-4501-4 Both probes feature two separate  $x_{short}$  and  $x_{long}$  functions to provide more accurate readings in the blue region over detectors

with only the  $x_{\mbox{\tiny long}}$  function and a simulated  $x_{short}$  using the z values. Sandwich filtering technology allows a very compact housing

design with a only 10 mm diameter diffuser window. This ensures a uniform illumination of the dif-



fuser as required for the measurement of illuminance.

A special feature of the CT-3701 detector is temperature control of it's photodiodes and filters for ambient temperature independent operation. This allows the

CT-3701 to be integrated into industrial qualification process control applications where long term stability is essential. CT-3701 and CT-4501 can be combined with other components like luminance steradian tubes.

Ordering In	formation &	Specificatio	ns											
Model	λresp	λresp	λresp	λresp	Т	ypical S	Sensitivi	ty	Sensing	cable	Opera-	Meter	plug	package
	Ch1	Ch2	Ch3	Ch4	CH1	CH2	CH3	CH4	Area	m	tion Temp.	Туре		page
	380-500 nm	500-730 nm	380-780 nm	380-580 nm		nA/lx at	2856 K	Ś	Diffuser		remp.			
CT-3701	x short	x long	у	Z	0.07	0.03	0.08	0.02	10 mm Ø	2	+40°C	P-9801	1	92
CT-4501	x short	x long	у	Z	0.07	0.03	0.04	0.01	10 mm Ø	2	5-40°C	X11	4	90
K-CT3701	Calibration c	f illuminance :	sensitivity A/lx,	x,y chromatio	city valı	ues and	color te	empera	ture. Incluc	les K-S	R in new o	detector (	order.	
K-CT4501	Calibration o	libration of illuminance sensitivity A/lx, x,y chromaticity values and color temperature. Includes K-SR in new detector order.												
KDW-S	Calibration o	libration of spectral sensitivity at one or multiple wavelengths without or in combination with accessory components												

# Luminance & Radiance

Luminance & Radiance

Radiance L is described as the intensity of optical radiation from an irradiated or irradiating surface in a particular direction measured with a specified solid angle.

The units of radiance are W/ (m2\*sr). Light detectors with a defined angle of view are used to measure radiance



**Photometric and Radiometric** 

detectors from the PD-16 series

combine with the LDM-9810 and

LDM-9811 viewing module to set

up a luminance or radiance

The PD-16VL01 photometric

detector has been adapted with

great care to the ideal V( $\lambda$ ) func-

tion, and in combination with the

I DM-98 meets the DIN class A specification for the adaptation

error of  $f_1$  ' less than or equal to

3%. Typical applications of the

luminance detector LDM-98/PD-

16VL01 include photometric

• Luminous surfaces such as

monitor screens, alphanumeric

displays, information signs,

emergency lighting and illumi-

• Reflecting surfaces such as

walls and workplace equip-

ment, projection screens, traf-

fic and information signs, road-

• Light emission in accordance

• Fluorescent pigments with

• In lighting design stage of

workplace and office design

with Federal Emissions Stat-

brightness measurement of:

nated control panels

ways and roller tracks

long afterglow decay

Luminance Measurement:

measurement device.

view.

Luminance L is described as the measurable photometric brightness of an illuminated or illuminating surface. In many cases luminance describes the quality of illumination much better than illuminance. The units of luminance are candela per unit area (cd/m<sup>2</sup>) or nit and footlamberts (fL). Detector heads

with a defined angle of view are used to measure luminance so that a small area within a larger uniform field is sampled.

> $cd/m^2 = nit = Im/(m^2 * sr)$  $fL = Im/(ft^2 * sr)$  $1 \text{ cd/m}^2 = 2.919 \text{x} 10^{-1} \text{ fl}$

#### LDM-9810 & LDM-9811: Luminance & Radiance Detectors

#### **Radiance Measurement:**

The PD-16RW05 radiometric detectors with spectral sensitivity within the wavelength range from 400 to 1000 nm are the right choice for radiance measurement in the visible and near infrared wavelength region.

#### **Light Hazard Protection:**

The quantity used for determination of blue light hazard and retinal thermal hazard to the human eye is radiance, measured in the radiometric unit of W/m<sup>2\*</sup>sr. Three different spectral response functions are typically used for light hazard protection radiance measurements:

- Blue-light Hazard (400-520 nm)
- Retinal Thermal Hazard (520-1400 nm)
- Retinal Thermal Hazard IR-A (780-1400 nm)

Gigahertz-Optik offers three different PD-16 detectors to be combined with the LDM-9811 viewer module for this application.

• PD-16BLH:

Blue light hazard detector

- PD-16RTH: Retinal thermal hazard detector
- PD-16RTHA: Retinal thermal hazard IR-A detector

#### **Optional Accessories:**

The LDM-98NL narrow lens provides a variable 1:1 to 2:1 near field magnification for measurement of spot diameters down to 0.1 mm.

The LDM-98FC enables coupling a flexible fibre bundle instead of a PD-16 detector to the LDM-98 to allow connection and use with spectral radiometers. The BHO-01 hard carrying case holds the LDM-98/PD16 with P-



9710 optometer for mobile applications.

#### Calibration in Light Measurement Units:

Calibration of the detector is carried out by Gigahertz-Optik's Calibration Laboratory for Optical Radiation Quantities, traceable to international standards.

Calibration in the units related to the selected PD-16 detector is supplied for each viewing angle with and without optional narrow lens.

For best handling of the multiple calibration factors the P-9710 and P-2000 optometers are recommended for use with the LDM-98/PD16 luminance and radiance detector heads outfitted with (-2) programmable data connectors



PD-16 Light Detector

LDM-9810 and LDM-9811 are modularly designed viewer modules that are combined with PD-16 detectors to form luminance and radiance measurement heads offering flexible spot size diameter selection. Both units are built with stable aluminium housings to ensure precise measurement in laboratory or field use.

The unit's wide aperture focussing objective is achromatically corrected, and satisfies the tightest requirements for the suppression of stray light and image formation.

Both viewer modules offer three selectable fields of view:

- LDM-9810: 20', 1° and 6°
- LDM-9811 .

1.7, 11 and 100 mrad

A f.o.v. selector knob is located on the rear of the device. The resultant size measurement spot is visible through the ocular viewfinder. In combination with the cross-hair targeting aid alignment of the LDM-98 to the zone of interest is simple.

Focussing is achieved by adjusting the sharpness of the image on the viewfinder screen.

The useable measurement distance of the LDM-98 spans from 0.3m to ∞. The area of the measured surface depends on the distance to the target and on the angle chosen for the field of

> Eve Eyepiece Surrounding Field Cross LDM-98 Hairs Detector Aperture Lens

utes

sample Beam Splitter Focusing Adjustment Measurement Spot

# LDM-9810 & LDM-9811: Luminance & Radiance Detectors









PD-16RW05

PD-16RTHA



PD-16VL01

F.O.V.		LDM-9810			LDM-9811	
Distance	<b>6</b> °	1°	20'	100 mrad	11 mrad	1.7 mrad
0.3 m	21	3.5	1.2	19.9	2.2	0.4
0.5 m	41.9	7	2.3	39.8	4.4	0.8
1 m	94.3	15.8	5.2	89.6	9.9	1.8
2 m	199	33.2	11	189.2	20.9	3.8
10 m	1037	173	57	986	109	19.8
th LDM-98NL Nar	row Lens (Optional)			L		
100 mm / f=∞	10.5	1.8	0.6	10	1.1	0.2
75 mm / f=0.3m	5.2	0.9	0.3	5	0.6	0.1

#### **Ordering Information & typical Specifications**

Model	λresp	Wavelength	f1	Тур	oical Sensitivi	ty	Imax	Sensing Area	cable		plug	packag
		Range	$\leq \%$	*) nA/cd*	*m² or **) nA	/W/m²sr	mA	Lens	m	Temp.		page
LDM-9810				20'	1°	6°		50 mm Ø		5-40°C		92
PD-16VL01	V(λ)	380-400	3	0.0035 *)	0.033 *)	0.9 *)	1	5.8 x 5.8	2	5-40°C	1,2,4	92
LDM-9811				1.7 mrad	11 mrad	100 mrad		50 mm Ø		5-40°C		92
PD-16BLH	BLH	400-520	-	0.2 **)	5.7 **)	380 **)	1	5.8 x 5.8	2	5-40°C	1,2,4	92
PD-16RTH	RTH	500-1200	-	0.3 **)	11 **)	780 **)	1	5.8 x 5.8	2	5-40°C	1,2,4	92
PD-16RTHA	RTHA	800-1200	-	0.2 **)	5 **)	360 **)	1	5.8 x 5.8	2	5-40°C	1,2,4	92
PD-16RW05	RW05	400-1000	-	0.38 **)	12 **)	840 **)	1	5.8 x 5.8	2	5-40°C	1,2,4	92
PD-16VL01	V(λ)	380-400	3	0.00004 *)	0.00015 *)	1.1 *)	1	5.8 x 5.8	2	5-40°C	1,2,4	92
LDM-98Z-NL	Narrow Lens	s for use with L	DM-9	810 and LDN	Л-9811							
LDM-98Z-FC	Fiber Couple	er to adapt flex	ible liq	ght guide bu	ndles to the L	DM-9810 a	and LD	M-9811				
K-LDM98BLH	Calibration o	of effective BLH	l radia	ance sensitiv	ity with <b>or</b> wit	hout LDM-9	98NL ir	n A/W/m²*sr (at	3 f.o.v.	angles). 1)		
K-LDM98RTH	Calibration of	of effective RTH	l radia	ance sensitiv	ity with <b>or</b> wit	hout LDM-	98NL ir	n A/W/m²*sr (at	3 viewir	ng angles). 1	)	
K-LDM98RTHA	Calibration of	of effective RT	HA ra	diance sensi	tivity with <b>or</b> w	without LDN	Л-98NL	_ in A/W/m²*sr (a	at 3 viev	ving angles).	. 1)	
K-LDM98RW05	Calibration of	of 400-1000 nm	n radia	ance sensitiv	ity with <b>or</b> wit	hout LDM-9	98NL ir	n A/W/m²*sr (at	3 viewir	ng angles). 1	)	
K-LDM98VL01	Calibration of	of luminance se	ensitiv	ity with <b>or</b> wi	thout LDM-98	BNL in A/cd	l/m² (a	t 3 viewing angle	es). 1)			
K-LDM98NL	Additional ca	alibration of lur	ninan	ce or radiand	ce sensitivity f	for LDM-98	with L	DM-98NL narro	w lens			
KDW-S	Calibration of	of spectral sens	sitivity	at one or mi	ultiple wavele	engths with	out or	in combination v	vith acc	essory comp	oonents	3
1) Includes K-SR	with new dete	ector order										

1) Includes K-SR with new detector order

Luminance

#### LDM-9901: Luminance Detector

Luminance detectors that are characterized as quality class A or B, in accordance with DIN 5032 Part 7, need to be constructed from carefully selected elements. These include optics with achromatic correction and low stray light, and a photopic detector with a response curve that has been carefully adapted to the ideal V( $\lambda$ ) function. Add an expensive telescopic finder and a focusing optical system for locating the target area and these luminance detector heads become too costly for many applications. The widely used auxiliary lenses that can be attached to illuminance detector heads for luminance measurement are not a good alternative, since they typically do not satisfy the requirements of DIN 5032 Part 7.

An economical alternative which does satisfy the requirements of IEC 61223-2-5 and DIN 5032 Part 7 is the LDM-9901. It is a high quality photometric detector for determining luminance in cd/ m<sup>2</sup>. Its properties satisfy the requirements for Quality Class B of DIN 5032 Part 7.

Distance t	to Spot Comparison
distance	spot diameter
0.5 m	~ 31 mm
0.7 m	~ 35 mm
1 m	~41 mm
5 m	~ 120 mm
10 m	~ 220 mm
50 m	~ 1000 mm
100 m	~ 2000 mm

# The LDM-9901 detector key components are:

- ©achromatically corrected low stray light lens to collect light with a 1° field of view
- Chigh-quality photodiode with accurate V( $\lambda$ ) correction filter
- ©sight for locating the target measuring area

The low-stray-light achromatic corrected objective is designed so that in its useable measurement range of 40 cm to infinity it is not necessary to make focusing adjustments to the lens. The 1° field of view means that the diameter of the measurement spot differs depending on the measuring distance.

The adaptation for photopic vision offers a f<sub>1</sub>'of  $\leq$  5 %. This means that the LDM-9901 detector can be used to measure the output of monitors, LED displays and other sources without suffering from an unacceptable increase in measurement uncertainty due to the quasimonochromatic type radiation of these sources.



троа	Adapter	(Optioi

LDM-9901 Comparison to DIN 5032 Cla	ass B Limi	ts (%)	
Characteristics	Symbol	LDM-9901	DIN
Calibration Uncertainty	U <sub>kal</sub>	1.5	4
$V(\lambda)$ Match Characteristic	f <sub>1</sub>	5	6
UV Response Characteristic	u	0.01	2
IR Response Characteristic	r	0.01	2
Linearity Characteristic	f <sub>3</sub>	0.2	2
Fatigue Characteristic (at 1 klx)	f <sub>5</sub>	0.1	1
Temperature Dependence Characteristic	f <sub>6</sub>	1	10



The LDM-9901 detector has an externally mounted targeting These back and front sights are adjusted for measuring distances of 70 cm, 140 cm and  $\infty$ as shown below. In order to locate the object to be measured, the detector head is held in front of the eye at such a distance (about 25 cm) that the back and front sights appear sharp. Depending on the observer's age, other distances may be appropriate. The two sights are then aligned with the desired part of the object being measured. The sights consists of a black

70 cm

slotted V-notch and a white ball. The following figure shows the different alignments for various measuring distances.

In many applications a luminance detector is used with a tripod. An adapter plate to attach the detector and tripod are offered as options.

The detector is connected to the optometer via a 1 m flexible connecting cable .

Calibration of luminance sensitivity in A/(m<sup>2</sup>cd<sup>-1</sup>) is carried out at Gigahertz-Optik's Calibration Laboratory and is confirmed by a factory certificate.

Infinity



140 cm

Notch and Bead Sight at difference Distances

Ordering Info	Drdering Information & typical Specifications												
Model	λresp	Wavelength f1		Typical Sensitivity	Imax	Sensing Area	cable	Operation	plug	package			
	Photometric	Range	≤%	pA/cd*m <sup>2</sup>	mA	Lens	m	Temp.		page			
LDM-9901	$V(\lambda)$	380-780 nm	5	25	1	22 mm Ø	2	0-40°C	1, 2, 4	92			
LDM-99Z-02	Ambient light	t shade made b	y elastic	rubber to place the LDM-99 d	lirect on	the monitor surfa	ce						
K-LDM9901	Calibration o	f integral lumina	ince sen	sitivity in A/cd/m². Including K	SR in n	new detector orde	r.						
KDW-S	Calibration o	f spectral sensit	ivity at or	ne or multiples wavelengths									

Ordering Information & typical Specification

# Luminance & Radiance / Luminous Flux & Radiant Power

#### SRT-M37: Front Lenses for 37-Type and 45-Type Illuminance & Irradiance Detectors



In some applications simple front lenses are required to limit the field of view of irradiance and luminance light detectors. The SRT-M37L series front lenses can be combined with the 37-type and 45-type irradiance,

illuminance and luminous color detectors for this purpose. The detector itself is mounted into the SRT-M45/37 base mount that has a M37x1 front thread where the SRT-M37L front lens can be attached.

The base mount offers a side M6 and camera type (1/4"-20 BSW) threaded hole for standard tripod or bench post or to fixture the device into an application. Beside the four front lenses offered custom solutions with quartz lenses, other lens diameters are available on special request.

Contact an applications engineer to discuss your requirements.

Ordering Information &	·· ·	1			1			1	
Model	λresp	Wavelength	Typica	al Sensitivity	Sensing Area	cable	Operation	plug	package
		Range			Lens	m	Temp.		page
SRT-M37L-1/VL-3701	$V(\lambda)$	380-780	1.2	pA/cd*m <sup>2</sup>		2	5-40°C	1, 2 or 4	89
SRT-M37L-2/VL-3701	V(λ)	380-780	2.3	pA/cd*m <sup>2</sup>		2	5-40°C	1, 2 or 4	89
SRT-M37L-5/VL-3701	$V(\lambda)$	380-780	12	pA/cd*m <sup>2</sup>		2	5-40°C	1, 2 or 4	89
SRT-M37L-1/CT-4501	V(λ)	380-780	0.05	pA/cd*m <sup>2</sup>		2	5-40°C	1, 2 or 4	89
SRT-M37L-2/CT-4501	V(λ)	380-780	0.4	pA/cd*m <sup>2</sup>	22.4 mm Ø	2	5-40°C	1, 2 or 4	89
SRT-M37L-5/CT-4501	V(λ)	380-780	2	pA/cd*m <sup>2</sup>	-	2	5-40°C	1, 2 or 4	89
SRT-M37L-10/CT-4501	V(λ)	380-780	9	pA/cd*m <sup>2</sup>	-	2	5-40°C	1, 2 or 4	89
SRT-M37L-1/RW-3703	Radiometric	400-800	5	nA/W/m²sr	-	2	5-40°C	1, 2 or 4	89
SRT-M37L-2/RW-3703	Radiometric	400-800	9	nA/W/m²sr	-	2	5-40°C	1, 2 or 4	89
SRT-M37L-5/RW-3703	Radiometric	400-800	48	nA/W/m²sr	-	2	5-40°C	1, 2 or 4	89
Specifications for other De	etector/Front Lens	Combinations	not shown	are available or	request. Quartz	lenses for	UV applicatio	ns are ava	ilable.
KDW-P3	Calibration of int	egral luminance	e sensitivity	of an illuminand	ce detector with f	ront lens ir	n A/cd/m²		
KDW-R3	Calibration of int	egral radiance s	sensitivity o	of an irradiance (	detector with fror	t lens in A/	W/m²*sr		
KDW-S3	Calibration of sp	ectral radiance	sensitivity	of an irradiance	detector with from	nt lens in A	/W/m²*sr nm		
SRT-M45/37B	Basic holder to r	nount 37-type a	nd 45-type	e detectors to SF	RT-M37 type com	ponents			
SRT-M37L-1	Front lens with 1	° field of view. N	leeds SRT	-M45/37B to cou	ple with detecto	r. 1 mm dia	ameter output	aperture.	
SRT-M37L-2	Front lens with 2	° field of view. N	leeds SRT	-M45/37B to cou	ple with detecto	r. 2 mm dia	ameter output	aperture.	
SRT-M37L-5	Front lens with 5	° field of view. N	leeds SRT	-M45/37B to cou	uple with detecto	r. 4.8 mm o	diameter outp	ut aperture	Э.
SRT-M37L-10	Front lens with 1	0° field of view.	Needs SR	T-M45/37B to co	ouple with detect	or. 10 mm	diameter outp	out apertur	e.
SRT-M37Z01	Ambient light sha	ade made by el	astic rubbe	er to be used wit	h the SRT-M37-L	lenses or	the SRT-M37-	-5 tube.	

#### **Luminous Flux & Radiant Power**



Luminous flux ♦ is the light power of any light source and is measured in lumens (Im). Since light is not typically emitted as a more or less parallel light beam, luminous flux must be measured with a measurement geometry that acquires the luminous flux independent of its spatial distribution. Light collecting integrating spheres or goniophotometers which can sum the angular light output over 360° are often used.

The luminous flux of a 100 W incandescent bulb is about 1380 lm, while that of a 20 W compact fluorescent lamp is about 1200 lm.

**Radiant power**  $\phi$  is the total power of optical radiation emit-

ted by a source measured in watts (W). The same measurement geometry as in luminous flux measurement applies only the detection system is not photopically corrected.

The most common method **to measure luminous flux and radiant power** from all directional emitting divergent beam sources employs an integrating sphere light detection system. Omnidirectional emitting light sources are placed in the center of the sphere, whereas the light from spot lamps enters the sphere through a port in the



Diffuse multiple reflections at the inner sphere wall suppresses disturbance effects created by the angle of inci-

dence, formation of shadows, reflections, modes etc.

#### sphere wall.

In the "Integrating Spheres" chapter several spheres from 8 mm up to 500 mm diameter are described which can be combined with Gigahertz-Optik light detectors and optometers to form complete luminous flux and radiant power light measurement systems.

# Luminous Flux & Radiant Power

#### LSM-9901: Luminous Flux Detector

tor is positioned close to the entrance port in accordance with DIN guidelines and allows for an acceptance angle of up to 90°. The detector is shielded from exposure to direct radiation by a baffle

The specifications of the photodiode correspond to DIN quality class B, permitting use with cold light sources, lamps with a high color temperature, or quasi-monochromatic radiation such as that from light emitting diodes.

The 50 mm diameter integrating sphere is coated with BaSO<sub>4</sub>, guaranteeing diffuse reflectivity and a spectrally neutral reflection characteristic. Optometer connections are made by means of a 2 m coaxial cable with appropriate plug type. Calibration of

luminous flux in lumens (Im) is carried out at Gigahertz-Optik's Calibration Laboratory for Optical Radiation Measurement, and is confirmed by a factory certificate.

Ordering Info	rmation & ty	pical Specifica	tions								
Model	λresp	Wavelength	f1	Typical Sensitivity	Imax	Sp	here	cable	1	plug	package
	Photometric	Range	$\leq$ %	nA/Im	mA	Port	Diameter	m	Temp.		page
LSM-9901	V(λ)	380-780 nm	5	40	1	12.5 mm Ø	50 mm	2	5-40°C	1, 2, 4	91
K-LSM9901	Calibration o	f integral luminc	us flux	sensitivity in A/Im. Ir	ncludes	K-SR in with I	new detector or	der.			
KDW-S	Calibration o	f spectral sensit	ivity at o	one or multiple wave	elength	s without or in	combination w	ith acce	ssory compo	onents.	



Integrating spheres are commonly used in radiant power measurement systems. Functioning as an optical integrator, the integrating sphere is able to measure radiation entering it's port independent of angle. The PRW-0505 radiometric detector for radiant power measurement in the visible and near infrared (NIR) is fitted with an integrating sphere. The sphere, with an internal diameter of 50 mm, is coated with Barium Sulfate offering a diffuse and spectrally flat reflectance characteristic within a specified wavelength range. BaSO4 offers a high re-

# PRW-0505: Polychromatic Radiant Power Detector

flectance of  $\sim$  97 % ensuring good sensitivity and also good long term stability.

The detector is positioned close to the sphere's entrance port to allow highly divergent radiation beams of up to 90 degrees to be accurately measured. The input port has a diameter of 12.7 mm. The radiometric detector is shielded from exposure to direct radiation from the measurement port by a baffle.

Its radiometric sensitivity covers wavelengths between 400-1000 nm allowing the radiant power of broadband and of quasi-monochromatic lamps (such as light emitting diodes) to be measured irregardless of their emission spectra. Gigahertz-Optik optometers are connected with a 2 m coaxial cable with an appropriate plug type. Calibration of sensitivity to radiant power



in W 400-1000 nm is carried out at Gigahertz-Optik's Calibration Laboratory for Optical Radiation Quantities, and is supported by a factory certificate.

Ordering Inf	Ordering Information & typical Specifications												
Model	λresp	Wavelength	Typical Sensitivity	фтах	Sph	ere	cable	Operation	plug	package			
	Radiometric	Range	μA/W	mW	Port	Diameter	m	Temp.		page			
PRW-0505	Polychromatic	400-1000 nm	90	500	12.5 mm Ø	50 mm	2	5-40°C	1, 2, 4	91			
K-LSM9901	Calibration of ir	ntegral radiant p	ower sensitivity in A	W. Includes K-S	R in new deteo	ctor order.							
KDW-S1	Calibration of s	pectral radiant p	oower sensitivity at c	one or multiple wa	avelengths wit	hout or in co	mbinati	on with acce	ssory co	mponents.			



**The LSM-9901** detector is used to measure the luminous flux, in lumens, of directional light beams with varying divergences. The measurement aperture size of its 50 mm diameter integrating sphere is 12.7 mm in diameter, which is sufficient for applications with endoscopes, light emitting diodes, optical fiber bundles and other beam emitters. The photometric light detec-

### LP-9901: Laser Power and Laser Stray Light Detector



W/m<sup>2</sup>

The LP-9901 Detector is designed to measure both laser power in W (the laser beam under fills the detector active area) and laser irradiance in mW/m<sup>2</sup> (the laser beam overfills the detector active area).

The measurement aperture has a diameter of 7 mm which correlates to the maximum size of the eve's pupil.

The LP-9901 offers a dynamic range from 0.1 µW to 100 mW (at 662.8 nm) for power measurement, and 0.0026  $\mu$ W/cm<sup>2</sup> to



260 mW/cm<sup>2</sup>, for the determination of maximum permissible exposure. The low profile detector can be securely held in the laser beam with a rigid 100 mm long detector handle built around the cable and connects to the optometer by a 2 m long cable

with appropriate plug type. Calibration of spectral sensitivity from 400 to 1100 nm is performed in 10 nm increments. Calibration is carried out at Gigahertz-Optik's Calibration Lab and is confirmed by a works certificate

Ordering In	formation & typ	ical Specificati	ons								
Model	λresp	Wavelength	Typical Se	ensitivity	<b>¢</b> ma	ах	cable	lmax.	Operation	plug	package
	Photodiode	Range	633 nm	900 nm	633 nm	900 nm	m	mA	Temp.		page
LP-9901	Si & ND Filter	400-1100 nm	1.3 mA/W	20 mA/W	100 mW	50 mW	2	1	0-40°C	1,2,4	91
K-LP9901	Calibration of sp	ectral radiant po	wer sensitivity	in A/W nm	and calculate	d spectral irr	adiance	sensitivit	y A/W/m²		
KDW-S1	Calibration of sp	ectral radiant po	wer sensitivity	at one or m	ultiple wavele	ngths in cor	nbinatio	n with acc	cessory comp	onents	

#### LP-01 & LP-02: Radiant Power Detectors with OP.DI.MA. Integrating Sphere





LP-01 detectors are designed for laser power and general optical power measurements in telecommunication testing systems. Using an integrating sphere



Typical spectral sensitivity LP-0101

size, high shunt resistance photodiodes providing high sensitivity and fast rise time to be used. As a result of multiple diffuse

reflectance, integrating spheres can reduce polarization effects, beam misalignment risk, signal bounce-back and PTD saturation. The machined OP.DI.MA. (optically diffuse material) spheres offer the highest reflectance (low attenuation) and longest term stability currently available. The LP-0101, LP-0102 and LP-0103 are built around a 30 mm diameter sphere with a 5 mm measurement aperture. A unique Gigahertz-Optik baffle design offers a large light acceptance angle with no risk of direct detector irradiation.

The LP-0201 employs an 8 mm diameter ODM integrating sphere with a 2 mm measurement port diameter for lowest attenuation. A low profile fast Si

<u>e</u>

photodiode is used enabling a very short rise time.

Common features include:

Diode array spectrometer can be coupled to an additional detector port with SMA-type fiber connector

Open port configuration is standard for direct measurement of laser diodes, LEDs or lasers.

Fiber connector adapters for FC, SC, ST and SMA connectors are available for the 30 mm sphere measurement port.

Calibration of spectral radiant power sensitivity in A/W nm is offered within the sensitivity range of each detector.



Typical spectral sensitivity LP-0103





**Laser Radiant Power** 

#### LP-01 Series: Radiant Power Detectors with OP.DI.MA. Integrating Sphere

Ordering Info	rmation & typi	cal Specificatio	ns									
Model	λresp	Wavelength	Typical Se	ensitivity	Spl	here	cable	Imax.	φmax	Operation	plug	package
	Photodiode	Range	633 nm	1300 nm	Port	Diameter	m	mA	mW	Temp.		page
LP-0101	InGaAs	800-1800 nm	-	0.4 mA/W	5 mm Ø	UPK-30-L	2	1	500	5-40°C	1,2,4	91
LP-0102	Si & InGaAs	350-1800 nm	0.6 A/W	0.4 mA/W	5 mm Ø	UPK-30-L	2	1	500	5-40°C	1,2,4	91
LP-0103	Si	350-1100 nm	0.6 mA/W	-	5 mm Ø	UPK-30-L	2	1	500	5-40°C	1,2,4	91
LP-0201	Si LP	400-1100 nm	30 mA/W	-	2 mm Ø	UPK-8-L	2	1	100	5-40°C	1,2,4	91
LP-01Z-01	Additional dete	ctor port for UP	<-30-L integra	ting sphere	of LP-01 for	r PD-11 type	detect	ors or U	FC-11/SI	MA socket		
UFC-11/FC	FC-type fiber s	ocket adapter fo	r the measure	ement port o	f LP-01 dete	ectors with L	JPK-30-	L integr	ating sph	nere		
UFC-11/SC	SC-type fiber s	ocket adapter fo	r the measure	ement port c	f LP-01 det	ectors with L	JPK-30	-L integr	ating spl	nere		
UFC-11/ST	ST-type fiber so	ocket adapter fo	r the measure	ment port o	f LP-01 dete	ectors with U	IPK-30-	L integr	ating sph	nere		
UFC-11/SMA	SMA-type fiber	socket adapter	for the measu	irement port	of LP-01 de	etectors with	UPK-3	0-L inte	grating s	phere		
K-LP01xx	Calibration of s	pectral radiant p	ower sensitivi	ity in A/W nn	n. xx = de	tector mode	l numb	er				
K-LP0201	Calibration of s	pectral radiant p	ower sensitivi	ity in A/W nn	n							
KDW-S1	Calibration of s	pectral radiant p	ower sensitivi	ity at one or	a chain of v	wavelength i	n comb	ination	with acce	essory		

#### LP-99 Series: Radiant Power Detectors with Barium Sulfate Integrating Sphere

**The LP-99** detectors are designed for general radiant power measurements of laser diodes and lasers in the wavelength range from 400 to 1700 nm.

The use of a 50 mm diameter integrating sphere as light collector offers a large 12.7 diameter measurement port and negates the need for using large diameter photodiodes. Small size photodiodes offer lower cost, higher shunt resistance and low capacitance. Also through its multiple reflectance characteristic, integrating spheres can reduce polarization effects, the risk of beam misalignment, signal bounce-back and PTD satura-



tion. Barium sulfate is a cost effective white diffuse coating for larger diameter integrating spheres. It's 97 % reflectance in the visible spectral range exhibits low attenuation making high sensitivity sphere based detectors

1,00

0,50

0,00

sensitivity

æ

are built around a 50 mm diameter sphere with a 12.5 mm aperture. A unique Gigahertz-Optik baffle design offers a large light acceptance angle without direct detector irradiation.

800 900 1000 1100 1200 1300 1400 1500 1600 1700 1800

wavelength(nm)

LP-9920



The LP-9910 is supplied with a Si photodiode covering the 400 to 1100 nm wavelength range.

The LP-9920 is supplied with an InGaAs photodiode for the 800 to 1800 nm range.

The LP-9930 offers both a Si and InGaAs photodiode for broadband use from 400 to 1800 nm with a 2-channel optometer.

The open port configuration is standard for direct measurement of laser diodes, LEDs or lasers. **Calibration** of spectral radiant power sensitivity in A/W nm is provided.



LP-9910

Ordering In	rdering Information & typical Specifications												
Model	λresp	Wavelength	Typical Se	ensitivity	Sph	ere	cable	lmax.	φmax	Operation	plug	package	
	Photodiode	Range	633 nm	1300 nm	Port	Diameter	m	mA	mW	Temp.		page	
LP-9910	Si	400-1100 nm	2.5 mA/W	-	12.5 mm Ø	UPB-50-L	2	1	500	5-40°C	1,2,4	91	
LP-9920	InGaAs	800-1800 nm	-	0.3 mA/W	12.5 mm Ø	UPB-50-L	2	1	500	5-40°C	1,2,4	91	
LP-9930	Si & InGaAs	400-1800 nm	2.5 mA/W	0.3 mA/W	12.5 mm Ø	UPB-50-L	2	1	500	5-40°C	1,2,4	91	
LP-99Z-01	Additional dete	ector port for UF	B-50-L integr	rating spher	e of LP-99 for	PD-11 type of	detecto	rs or UF	-C-11/SN	1A socket		-	
K-LP99xx	Calibration of s	spectral radiant	power sensiti	ivity in A/W r	nm. xx = det	ector model	numbe	r					
KDW-S1	Calibration of s	spectral radiant	power sensiti	ivity at one o	or multiple way	elengths in o	combir	ation wi	th acces	sory compor	nents		

# Laser Power & Laser Pulse Shape

#### LPPA-99: Laser Power and Pulse Shape Detector



Optoelectronic sensing systems such as laser range finders, profile scanners, distance sensors for example, make use of pulsed lasers generating high peak power and pulse lengths in the order of nanoseconds.

For general analysis or for reasons of safety the pulse shape of a single pulse or pulse train needs to be measured, a high speed detector with a rise time shorter than that of the laser pulse is required.

These photodiodes must be very low capacitance as well as offer

Tek Run: 2.00GS/s Average

C1 Rise 4.385ns C1 High 3.14mV C1 Rise 4.385ns C1 High 3.14mV H 10.0ns Ch1 7 1.22mV 19 Apr 2001 11:00:32

a deep depletion zone to avoid photon absorption outside the depletion zone.

Photodiodes with small active areas cannot be used to measure radiant power at a risk that the laser beam would be larger than the photodiodes area creating measurement error.

Conventional PIN photodiodes require a large bias voltage for full depletion for long wavelengths above 800 nm.

The LPPA-99 series detectors are designed around a 50 mm diameter integrating sphere. The



large 12.7 mm diameter measurement port allows the measurement of large diameter laser beams and laser diodes with divergent beams. A low-profile type Si photodiode provides a fast rise time even at 1060 nm. A baffle protects the photodiode from direct radiation exposure from the entrance port.

These unique features make the LPPA-99 detectors the right choice for pulse analysis of short and high peak power laser pulses in the absolute units of watts (W).

The Low Profile photodiode is operated with 12 V reverse bias, and is placed in series with a load resistor to create a voltage signal that correlates to the radiant power. The detector connects to a fast sensitive oscilloscopes to display the pulse shape.

Digital storage oscilloscopes enable the user to store individual pulse shapes.

Three LPPA-99 models are available with different sensitivities and rise times.

The spectral sensitivity of the detectors in V/W nm is calibrated in 10 nm increments from 400 to 1060 nm.

Calibration is carried out at Gigahertz-Optik's Calibration Laboratory for Optical Radiation Quantities, and is confirmed by a factory certificate.



Output: SMC Socket for 50 Ohm Coaxial Cable

Ordering In	Ordering Information & typical Specifications													
Model	λresp	Wavelength	Typical S	Sensitivity	<b>þ</b> min	<b>¢</b> max	Rise	Sphere Din	here Dimension		Load	plug	package	
	Photodiode	Range nm	633 nm	900 nm	@80	)0 nm	Time ns	Port	Diameter	m	Ohm		page	
LPPA-9901	Si LP	400-1060	50 mV/W	50 mV/W	5 mW	100 W 1)	≤ 8	12.5 mm Ø	50 mm	2	50	BNC	91	
LPPA-9902	Si LP	400-1060	15 mV/W	15 mV/W	15 mW	333 W 1)	≤5	12.5 mm Ø	50 mm	2	50	BNC	91	
LPPA-9903	Si LP	400-1060	1 mV/W	1 mV/W	200 mW	5 kW 1)	≤2	12.5 mm Ø	50 mm	2	50	BNC	91	
K-LPPA99 Calibration of spectral radiant power sensitivity in V/W nm.														
1) Pulse pea	k power only!	Average pov	ver should	not exceed	1 500 mW									

#### **Modular Light Detectors**

Modular light detectors offer a mechanical interface which allows accessory components such as optical filters, integrating spheres, diffuser windows and other front-end optics to be easily attached and interchanged. A wide range of base detectors are stocked of different photodiode type and active area. This permits selection of the best fit photodiode spectral sensitivity, light responsivity, rise and fall time for the specific job.

The following spectral plots show the typical spectral sensitivities of the photodiodes available in the PD-11, TD-11 and MD-37 series modular detectors. Note that spectral sensitivity will

change when the detector is combined with an input optic component including integrating spheres. For this reason the actual measured spectral sensitivity of any detector accessory combination is provided in the detector calibration certificate. Calibrations are optional and ordered separately.



The **PD-11** series modular detectors are primarily designed to be used in combination with Gigahertz-Optik integrating spheres. A typical application example for sphere based detectors is the measurement of radiant power of lasers, laser diodes and other monochromatic light sources.

**The PD-11** detectors offer a compact package with flange mount to allow direct mounting to all spheres with 11-type port adapters. Arrangement of discrete linear arrays with a grid size

of 11 mm are also possible. Calibrations are offered with or without accessory.



Ordering Information & typical Specifications													
Model			Photodiode	Э		Ту	o. Sensitivity	A/W	lmax.	Operating	plug	cable	package
	λresp	Resp.	$\lambda$ range nm	mm <sup>2</sup>	Size mm	@ 350 nm	@ 900 nm	@ 1350 nm	mA	Temp °C		m	page
PD-1101	Si	Page	200-1100	13	3.6 x 3.6	0.15	0,5	-	1	5 to 40	1,2,4	2	90
PD-1102	InGaAs	85	800-1800	0.07	0.3 Ø	-	0.14	0.9	1	5 to 40	1,2,4	2	90
PD-1103	InGaAs	85	800-1800	0.79	1Ø	-	0.14	0.9	1	5 to 40	1,2,4	2	90
PD-1104	InGaAs	85	800-1800	7.1	3 Ø	-	0.14	0.9	1	5 to 40	1,2,4	2	90
PD-1105	SiC	85	200-400	1	1.1 x 1.1	0.02	-	-	0.05	5 to 40	1,2,4	2	90
PD-1106	Si LP	85	400-1050	0.25	0.5 x 0.5	0.06	0.4	-	1	5 to 40	1,2,4	2	90
PD-1107	Si LP	85	400-1050	3.6	1.9 x 1.9	0.06	0.4	-	1	5 to 40	1,2,4	2	90
PD-1108	Si LP	85	400-1050	11.9	3.4 x 3.4	0.06	0.4	-	1	5 to 40	1,2,4	2	90
PD-1109	GaP	85	250-550	1.0	1.1 x 1.1	0.07	-	-	0.1	5 to 40	1,2,4	2	90
PD-1110	GaAsP	85	200-760	5.2	2.3 x 2.3	0.04	-	-	0.1	5 to 40	1,2,4	2	90
K-PD11xx	Calibration of relative spectral sensitivity. xx = detector model number												
KDW-S1	Calibration of spectral radiant power sensitivity at one or multiple wavelengths in combination with an integrating sphere												
KDW-S2	Calibration	Calibration of spectral irradiance sensitivity at one or multiple wavelengths in combination with a diffuser											
KDW-S3	Calibration	n of spec	ctral radiance	sensitivit	y at one or r	nultiple wav	elengths in co	ombination wi	th radia	nce front tul	ce/lens		
KDW-S4	Calibration	n of spec	ctral radiant ir	ntensity se	ensitivity at c	one or multip	le wavelengt	hs					

# Modular

# TD-1101, TD-1102 and TD-1103: Temperature Controlled Modular Detectors



The **TD-11** family of detectors are designed with a unique temperature control feature. The photodiode is precisely heated and maintained at a higher operation temperature than the ambient temperature. This function ensures measure-

ment accuracy when these detectors are to be used in applications with widely varying ambient temperatures as in laser power measurements involving sphere based detectors for example. The control electronics are housed in an external box with a plug in power supply.

Both silicon and indium gallium arsenide photodiodes are offered for the wavelength range

200 to 1100 nm and 800 to 1800 nm respectively.

TD-11 detectors can be direct mounted to all spheres with 11-type port adapters.

Calibrations are offered for the detector with or without any accessory components.

Ordering I	Ordering Information & Specifications												
Model			Photodiode	;		Тур	o. Sensitivity	A/W	Imax.	Operating	plug	cable	package
	λresp	page	$\lambda$ range nm	mm <sup>2</sup>	Size mm	@ 350 nm	@ 900 nm	@ 1350 nm	mA	Temp °C		m	page
TD-1101	Si	27	200-1100	13	3.6 x 3.6	0.15	0.5	-	1	55	1,2,4	2	90
TD-1102	InGaAs	27	800-1800	0.07	0.3 Ø	-	0.14	0.9	1	55	1,2,4	2	90
TD-1103	InGaAs	27	800-1800	0.79	1 Ø	-	0.14	0.9	1	55	1,2,4	2	90
K-TD11xx	Calibration of relative spectral sensitivity. xx = detector model number												
KDW-S1	Calibration of spectral radiant power sensitivity at one or multiple wavelengths in combination with an integrating sphere												
KDW-S2	Calibration of spectral irradiance sensitivity at one or multiple wavelengths with a diffuser												
KDW-S3	Calibration of spectral radiance sensitivity at one or multiple wavelengths in combination with radiance front tube/lens												
KDW-S4	Calibratior	n of spe	ectral radiant	intensity	sensitivity	at one or m	ultiple wav	elengths wav	relength	าร			

#### **MD-37: Modular Detector**



#### The modular light detectors of the MD-37 series feature a M30x1 threaded front interface for fast and simple screw-on assembly of mating accessory components.

For multi-purpose use modular detectors offer a wide range of base photodiode detector types and different active areas. This provides a great degree of freedom in selecting a spectral sensitivity, light responsivity, rise and fall time for a specific job. M30x1 is the standard thread size for all 37-type accessories which mate to the MD-37 detectors.

Some typical detector set-ups for the MD-37 series include:

 Broadband irradiance detectors (MD-37 with MD-37Z- DR11 diffuser front cap)

 Narrow band irradiance detectors (MD-37 with narrow band pass filter assembled in MD-37Z-H7 plus diffuser front cap MD-37Z-FR11)

Calibrations are offered with or without accessory components attached.

Ordering Inform	mation &	Specifi	cations										
Model			Photodiode	;		Ty	o. Sensitivity	A/W	Imax.	Operating	plug	cable	package
	λresp	Resp.	$\lambda$ range nm	mm <sup>2</sup>	Size mm	@ 350 nm	@ 900 nm	@ 1350 nm	mA	Temp °C		m	page
MD-37-GP6	GaP	27	250-550	6.25	2.5 x 2.5	0.07	-	-	0.1	5 to 40	1, 2, 4	2	90
MD-37-GAP5	GaAsP	27	200-760	5.29	2.3 x 2.3	0.04	-	-	0.1	5 to 40	1, 2, 4	2	90
MD-37-IGA1	InGaAs	27	800-1800	0.79	1 Ø	-	0.14	0.9	1	5 to 40	1, 2, 4	2	90
MD-37-IGA3	InGaAs	27	800-1800	7.1	3 Ø	-	0.14	0.9	1	5 to 40	1, 2, 4	2	90
MD-37-SC1	SiC	27	800-1800	1	1 sq	0.02	-	-	0.05	5 to 40	1, 2, 4	2	90
MD-37-SU13	Si	27	200-1100	12.96	3.6 x 3.6	0.15	0.5	-	1	5 to 40	1, 2, 4	2	90
MD-37-SU33	Si	27	200-1100	33,64	5.8 x 5.8	0.15	0.5	-	1	5 to 40	1, 2, 4	2	90
MD-37-SU100	Si	27	200-1100	100	10 x 10	0.15	0.5	-	1	5 to 40	1, 2, 4	2	90
MD-37Z-C11	Front cap	with 11	mm diamet	er aperti	ure								
MD-37Z-DR11	Front cap	o with 11	mm diamet	er RADI	N diffuser								
MD-37Z-H7	Front ada	apter an	d/or filter hol	der for C	).5 to 1 inch	n diameter fil	ters						
MD-37Z-H7xx	Front ada	apter wit	h custom sp	ecified b	band-pass	or edge filter	. xx = filter ty	rpe					
K-PD37xx	Calibratio	on of rela	ative spectra	l sensitiv	/ity.xx = d	letector mod	el number						
KDW-S1	Calibratio	on of spe	ectral radiant	powers	sensitivity a	at one or mul	tiple waveler	ngths in comb	inatior	n with an inte	egrating	sphere	;
KDW-S2	Calibratio	on of spe	ectral irradiar	nce sens	sitivity at on	ne or multiple	wavelength	s in combinat	ion wit	h a diffuser			
KDW-S3	Calibratio	on of spe	ectral radian	ce sensi	tivity at one	e or multiple	wavelengths	in combinatio	on with	radiance fr	ont tube,	/lens	
KDW-S4	Calibratio	on of spe	ectral radiant	intensit	y sensitivity	/ at one or m	ultiple wavel	engths					

🖲 Gigahertz-Optik

# Light Detectors - Page 87

Modular

#### PD-93: 1 cm<sup>2</sup> Circular Area Modular Light Detectors

Modular light detectors using photodiodes with a one square centimetre circular active area have been offered by Gigahertz-Optik since 1986. The symmetry of a circular detection area combined with a large 11.28 mm diameter has made these detectors a standard for many applications.

The area around the light sensitive 1 cm<sup>2</sup> area is metal shielded and thus insensitive to light so that irradiance (W/cm<sup>2</sup>) can eas-

1.00

0.00

measurements.

360

sensitivity 0'20

ē

ily be calculated from a known radiant power (W) assuming the active area is to be evenly overfilled.

The photodiodes mounted into a housing that can be fastened into optical test arrangements using the threaded M6 and 1/4-20 side holes. The front of the housing is internally threaded into which accessories can be attached (common 1.3/16-24-UNS thread size). In addition to two different photodiodes, Giga-

hertz-Optik offers accessories to form a wide variety of detector configurations The simple design of the mechanical parts also means that special custom fabrications can be prepared easily and economically for soapplications. Our sales engineers are ready to work with you and provide advice on standard and custom configurations.



PD-9304: 1 cm<sup>2</sup> Si-Photodiode

660 760 860

wavelength(nm)

The PD-9304 is a large, 1 cm<sup>2</sup> photodiode for

general light measurement purposes and for

use in combination with PD-93VL and PD-

93RW photometric and radiometric correction

filters for illuminance and integral irradiance

960 1060



**PD-93VL: Photometric Filter** 

Typical spectral sensitivity with PD-9304

The **PD-93VL** V( $\lambda$ ) correction filter screws directly onto the PD-9304 detector. It's built into the PD-9302 filter holder which is internally and externally threaded (1.3/16-24-UNS) for mating to PD detectors and to accept accessories.

**PD-93RW: Radiometric Filter** 



#### Typical spectral sensitivity with PD-9304

The **PD-93RW** filter built into the PD-9302 filter holder, gives the PD-9304 a radiometric flat spectral sensitivity within the wavelength range from 400 and 1000 nm.

13mm

6mm

#### PD-93COS: Diffuser

460 560

The **PD-93COS** diffuser cap generates a field of view approximating a cosine-corrected function, enabling the measurement of illuminance or irradiance. The diffuser screws onto the PD-9302 or PD-9304.

#### PD-9302: Universal

The **PD-9302** universal holder is 1.3/16-24-UNS internally and externally threaded suitable for attachment to the **PD-9301** detector housing. One inch diameter filters like neutral density, band pass and correction filters can be mounted into it.

#### PD-9303: Front Cap

The PD-9303 cap can be used as an end piece on the PD-9301

One of the parameters most commonly used to specify LEDs is the directional output intensity. This radiometric or photometric intensity is derived from the flux incident on a detector by the LED at specified distance.

A detector with a circular entrance aperture having 100 mm<sup>2</sup> detector housing, or the PD-9302 universal holder.

PD-9302: Modular Schematic

27.2mm

Photodiode PD-9301 with BNC Connector PD-9302 PD-9303

#### PD-9311 & PD-9312: Averaged LED Intensity Front Tubes

(11.3 mm diameter) is specified for the CIE standard A and B measurement conditions for average LED intensity. The distance between the front tip of the LED and the plane of the entrance aperture of the detector is stated at d = 316 mm for CIE standard condition A and d =



100 mm for standard condition B. This corresponds to solid angles of view of 0.001 sr for Condition A and 0.01 sr for Condition B. The PD-9311 and PD-9312 front tubes attach to the PD-9304 detector with either PD-93VL photometric or PD-93RW radiometric filter for CIE Condition A and B conformal

25mm

measurements. Calibration of the average LED Intensity sensitivity in cd and W/sr is offered by Gigahertz-Optics calibration laboratory.

# Cigahertz-Optik

# Modular, Weather Proof

# PD-9306: 1 cm<sup>2</sup> Circular Area Light Detector

The **PD-9306** is a large area 1 cm<sup>2</sup> UV-enhanced silicon photodiode for spectral radiant power and/or spectral irradiance measurements. It's large diameter active area offers higher signals for applications with low intensity monochromatic radiation.

The detector is supplied in the PD-9301 housing with internal 1.3/16-24-UNS front thread to accept other accessory compo-

nents. A 2 m long coaxial cable with BNC plugs can be coupled to the BNC connector at the end of the photodiode. Calibration of spectral radiant power sensitivity is supplied by the Gigahertz-Optik's calibration laboratory as required.



### **Ordering Information & Specifications**

Model	λresp	Wavelength	Typ	ical Sensitivity	Sensing Area	Imax	cable	Operating	plug	package			
Nodel	Niesp	Range		ical Sensitivity	Photodiode	mA	m	Temp °C	piug				
		0				IIIA	111			page			
PD-9304	Si Photodiode	380-1100 nm	0.4	A/W @ 633 nm	1 cm <sup>2</sup>	1	2	5 to 40	1,2,4	91			
PD-93VL&PD-9304	$V(\lambda)$	380-780 nm	28	nA/lx	1 cm <sup>2</sup>	1	2	5 to 40	1,2,4	91			
PD-93COS&PD-93VL&PD9304	$V(\lambda)$	380-780 nm	2.8	nA/lx	1 cm <sup>2</sup>	1	2	5 to 40	1,2,4	91			
PD-9311&PD93-VL&PD-9304	$V(\lambda)$	380-780 nm	2.6	μA/cd	1 cm <sup>2</sup>	1	2	5 to 40	1,2,4	91			
PD-93RW&PD-9304	Radiometric	400-1100 nm	0.2	A/W @ 800 nm	1 cm <sup>2</sup>	1	2	5 to 40	1,2,4	91			
PD-93COS&PD-93RW&PD-9304	Radiometric	400-1100 nm	16	mA/W/cm <sup>2</sup>	1 cm <sup>2</sup>	1	2	5 to 40	1,2,4	91			
PD-9311&PD-93RW&PD-9304	Radiometric	400-1100 nm	1.8	mA/W/sr	1 cm <sup>2</sup>	1	2	5 to 40	1,2,4	91			
PD-9306	UV Si	250-1100 nm	0.16	A/W @ 350 nm	1 cm <sup>2</sup>	0.2	2	5 to 40	1,2,4	91			
K-PD9304	Calibration of s	pectral radiant	power s	ensitivity in A/W/nr	n from 350 to 1 <sup>-</sup>	100 n	m in 10	) nm steps					
K-PD9306	Calibration of s	pectral radiant	power s	ensitivity in A/W/nr	n from 250 to 1 <sup>-</sup>	100 n	m in 10	) nm steps					
KDW-P2	Calibration of il	luminance sen	sitivity in	lx of PD-9304/PD-	93VL with or wit	hout	PD-930	COS					
KDW-P4	Calibration of a	Calibration of average LED Intensity sensitivity in cd of PD-9						D-9304&PD-93VL&PD-9311 or PD-9312					
KDW-R2	Calibration of irradiance sensitivity in W/m <sup>2</sup> of PD-9304/PD-93RW with or without PD-93COS												
KDW-R4	Calibration of average LED Intensity sensitivity in W/sr of PD-9304&PD-93RW&PD-9311 or PD-9312												
KDW-S	Calibration of s	pectral sensitiv	rity at on	e or multiple wave	elengths without	or wi	th acce	essory comp	onents	S			

#### **WPD: Weather Proof Detector**



The **WPD** series weather proof detectors are designed for year round outdoor light measurements. To withstand all environmental conditions the light sensor is shielded in a rugged metal housing. The detector aperture is enclosed and sealed under a precision made quartz dome. Photodiode, filter and diffuser are

temperature stabilized to +25 deg C by a pettier cooler. Temperature is stabilized to +/-1 deg C for operation independent of the ambient conditions. The specified operating temperature range is -20 to +50 deg. C. The electronics of the **PTC-0101** pettier temperature controller are supplied in a separate compact

housing for bench-top or rackmount use with cable connections.

In principle the sensor elements of all Gigahertz-Optik light detectors can be integrated into the WPD housing. The most common applications to date are in the ultraviolet spectral range. For example WPD detectors can serve as reference detectors for high resolution spectral radiometers monitoring global solar UV radiation The WPD water proof light detectors are normally connected with a Gigahertz-Optik signal amplifier or optometer for direct measurement of current or irradiance as well as data logging.

Due of the custom nature of this detector and supporting electronics interested customers should contact the factory to discuss his specific requirements.





**Light Detector Accessories** 

#### SRT-M37: Front Tubes and Front Lenses

**SRT-M45/37B:** Base holder for coupling 45-type and 37-type standard light detector and 37-type solid angle tubes and SRT-37L front lenses.

M30x1 front thread. M6 and camera tripod (1/2"-20BSW) side mounting holes.

#### Front Tubes:

**SRT-M37-50:** 37 mm outside diameter tube, 50 mm length internally threaded over its complete length and externally threaded front end (both M30 x 1). Can be machined to any length down to 12.5 mm.

SRT-37-AB Apertures can be glued into the tube at either end. SRT-37/25: Reducer to couple

SRT-M37-50 to SRT-M25-25 to decrease the outside diameter of solid angle tubes.

**SRT-25-25:** 25 mm outside diameter tube, 25 mm length internally threaded over its complete length and externally threaded front end (both M20 x1). Can be machined to any length down to 12.5 mm.

**SRT-37-AB:** Aperture plate insert for SRT-M37 tubes. 2 mm diameter centering hole for drilling the aperture opening to the desired diameter.

#### Front Lenses:

The SRT-M37L Front Lenses are used with the 37-type (eg. VL-3701, RW-3701) or 45-type (eg. HCT-99) Gigahertz-Optik detectors to form economical luminance and radiance measurement set-ups. The detector plus front lens can be used for relative measurement or can be calibrated for direct readout in cd/m<sup>2</sup> or W/m<sup>2</sup>/sr with a Gigahertz-Optik optometer.

The SRT-M37L set-up consists of the SRT-



M45/37B base holder which mounts onto the detector by means of its M30x1 threaded front and one of the SRT-M37L lenses mounted into it. Lenses are available with different fieldof-views. The 25 mm diameter BK7 lens provides a "flat" wave-

length range from 380 nm to 1500 nm. Quartz lenses are available on request



Ordering Information & Specifications								
Model	Description	Min. Required Detector Aperture	Specifications with Detectors on Page					
SRT-M45/37B	Base holder to couple 45 and 37 type light detectors to SRT-M37 accessory	-	-					
SRT-M37-50	37 mm diameter front tube, 50 mm length and M30x1 thread	-	-					
SRT-M37/M20	Adapter to mount 37 and 25 tubes together	-	-					
SRT-M20-25	25 mm diameter front tube, 25 mm length and M20x1 thread	-	-					
SRT-M37-AP	Aperture plate to glue into SRT-M37-50 front tube. 2 mm center hole	-	-					
SRT-M37L-1	Luminance & radiance front lens for SRT-M45/37B with 1° field-of-view	0.7 mm Ø	80					
SRT-M37L-2	Luminance & radiance front lens for SRT-M45/37B with 2° field-of-view	1.4 mm Ø	80					
SRT-M37L-5	Luminance & radiance front lens for SRT-M45/37B with 4.8° field-of-view	3.4 mm Ø	80					
SRT-M37L10	Luminance & radiance front lens for SRT-M45/37B with 10° field-of-view	7 mm Ø	80					



# Light Detector Accessories / Dimensions

/WG

/WQ

### /WQ: Waterproof Detector Option

	ing: © Oceanog © Limnolog © Aquacultu © Solar Biod	y	and other user specific require- ments are available through Gi- gahertz-Optik's custom design services. A transparent dome, made from glass or UV transmitting quartz covers the cosine diffuser. The detector housing is sealed and	filled with a drying agent (desiccant). This option cannot be retrofit as an after sales modi- fication! Please contact our application engineering department to dis- cuss your underwater light meas- urement project.
/WQ Waterproof Option e	Ordering Info	ormation & Specificat	tions	
ively seals any type-37 series	Model	Description		

nal modifications. Ordering P/N example: PS-3703-2/WG.

modifications. Ordering P/N example: PS-3702-2/WQ.

The fectively seals any type-37 s detector housing. This modification enables full submersion of the detector in water and other liquids. UV/Vis/NIR underwater detectors can be custom tailored

#### /LC: Cable Extension

For long distance applications requiring the standard 2 m cable between optometer and light

detector be extended. This option is offered for single channel detectors with BNC type (-1) and

calibration type (-2) detectors only. For a VL-3701 detector requiring a total cable length of 5

Waterproof modification of 37 type detectors. Glass dome cap for VISIBLE to NIR , and inter-

Waterproof modification of 37 type detectors. Quartz dome cap for UV to NIR, and internal

m three 3 m /LC should be ordered since 2 m is already included in the detector price.

# **Ordering Information & Specifications**

Model Description /LC

Cable extension in 1 m steps for light detectors with BNC type (-1) or calibration data type (-2) connectors

#### **BHO: Carrying Cases**

Carrying cases are offered for storage and transporting light

detectors and optometers. Please see Optometer chapter for information on available models



**11 Type Housing Dimensions** 

#### **37 Type Housing Dimensions**



C

**TD-11 Type Housing Dimensions** 



#### **45 Type Housing Dimensions**



### **MD-37 Type Housing Dimensions**





### RCH-0, RCH-1 and RCH-5 Housing Dimensions



LP-99, LSM-99 and PRW-05 Housing Dimensions



**LP-02 Housing Dimensions** 



Ø 15 mm Ø 15 mm Ø 15 mm Ø 15 mm Ø 16 mm Ø 63 mm Ø 6 mm

# **LP-01 Series Housing Dimensions**



### LPPA-9901 Housing Dimensions



# Page 92 - Light Detectors

# Light Detector Dimensions

